

**EE 472 Lab 4**  
**Learning the Development Environment - The Next Next Step**

Jonathan Ellington  
Patrick Ma  
Jarrett Gaddy

## Contents

1	Abstract	1
2	Introduction	1
3	Discussion of the Lab	1
3.1	Design Specification . . . . .	1
3.1.1	Specification Overview . . . . .	1
3.1.2	Identified Use Cases . . . . .	2
3.1.3	Detailed Specifications . . . . .	5
3.1.4	Detailed Task Specifications . . . . .	6
3.2	Software Implementation . . . . .	11
3.2.1	Functional Decomposition . . . . .	11
3.2.2	System Architecture . . . . .	12
3.2.3	High-level Implementation in C . . . . .	14
3.2.4	Task Implementation . . . . .	15
4	Presentation, Discussion, and Analysis of the Results	18
4.1	Results . . . . .	18
4.2	Discussion of Results . . . . .	19
4.3	Analysis of Any Errors . . . . .	19
4.4	Analysis of Implementation Issues and Workarounds . . . . .	20
5	Test Plan	21
5.1	Test Specification . . . . .	22
5.1.1	Scheduler . . . . .	22
5.1.2	Measure Task . . . . .	22
5.1.3	Compute Task . . . . .	23
5.1.4	Keypad Task . . . . .	23
5.1.5	Display Task . . . . .	23
5.1.6	Warning/Alarm Task . . . . .	23
5.1.7	Serial Task . . . . .	24
5.1.8	Status Task . . . . .	24
5.1.9	EKG Capture . . . . .	24
5.1.10	EKG Process . . . . .	24
5.1.11	Command Task . . . . .	24
5.1.12	Remote Communications . . . . .	25
5.2	Test Cases . . . . .	25
6	Summary and Conclusion	26
6.1	Final Summary . . . . .	26
6.2	Project Conclusions . . . . .	26
A	Breakdown of Lab Person-hours (Estimated)	28
B	Supported Web Browser Commands	29
C	Activity Diagrams	30

D	Source Code	37
D.1	Main Function . . . . .	37
D.2	Global Data . . . . .	46
D.3	Timebase . . . . .	49
D.4	Modified lmi_fs.c . . . . .	50
D.5	Circular Buffer . . . . .	56
D.6	Warm up File (For Interrupt Initialization) . . . . .	58
D.7	Tasks . . . . .	63
D.7.1	Task Control Blocks . . . . .	63
D.7.2	Startup Task . . . . .	64
D.7.3	Measure Task . . . . .	70
D.7.4	Compute Task . . . . .	75
D.7.5	Keypad Task . . . . .	77
D.7.6	Display Task . . . . .	81
D.7.7	Warning/Alarm Task . . . . .	85
D.7.8	Serial Task . . . . .	92
D.7.9	EKG Capture . . . . .	95
D.7.10	EKG Process . . . . .	98
D.7.11	Command Task . . . . .	101
D.7.12	Status . . . . .	111
D.8	Website Source . . . . .	112

## List of Tables

1	Specifications for measurement data . . . . .	6
2	Empirically determined task runtimes . . . . .	19
3	Estimated instructions per task, rounded to the nearest instruction . . . . .	19
4	Initial values and warning/alarm states . . . . .	24

## List of Figures

1	Use case diagram . . . . .	3
2	Expected screen layout of remote serial terminal . . . . .	7
3	Expected OLED display screens . . . . .	9
4	Expected Web Browser Interface . . . . .	11
5	Top level system functionality . . . . .	12
6	Low Level Functional Decomposition . . . . .	12
7	System Architecture Diagram . . . . .	14
8	Control Flow Diagram . . . . .	16
9	Measure Activity Diagram . . . . .	30
10	Compute Activity Diagram . . . . .	31
11	Keypad Activity Diagram . . . . .	32
12	Warning Activity Diagram . . . . .	33
13	Serial Activity Diagram . . . . .	33
14	Display Activity Diagram . . . . .	34
15	Status Activity Diagram . . . . .	34
16	EKG Capture Activity Diagram . . . . .	35
17	EKG Process Activity Diagram . . . . .	35

18	Command Activity Diagram . . . . .	35
19	Remote Activity Diagram . . . . .	36

## 1 ABSTRACT

In this lab the students are to take on the role of an embedded system design team. They will design modifications to the medical instrument previously designed. When the device finds metrics are out of the acceptable range, the user is notified, thus saving them from potential health risks. The students first laid out the design for their system using various design tools, then they implemented the system in software. Finally the students tested their system and verified that it is ready to start saving lives.

## 2 INTRODUCTION

The students are to design an embedded system on the Texas Instruments Stellaris EKI-LM3S8962 and EE 472 embedded design testboard. The design must implement a medical monitoring device. This device must monitor a patient's temperature, heart rate, and blood pressure, as well as its own battery state. The design must indicate when a monitored value is outside of a specified range by flashing an LED on the test board. When a value deviates even further from the valid range an alarm will sound. This alarm will sound until the values return to the valid range or the user acknowledges the alarm with a button. The values of each measurement will also be printed to the OLED screen. This implementation will build upon the previous implementation of the device by adding functionality. Added functions include heart rate sensor, keypad input, menu display, data logging, and UART serial communication.

The design will be tested to verify proper behavior on alarm and warning notifications. In addition the implementation will be tested by measuring the amount of time that each of the 8 program tasks running the instrument take to execute. These tasks are mini programs that each handle a part of the instruments purpose.

## 3 DISCUSSION OF THE LAB

### 3.1 Design Specification

#### 3.1.1 *Specification Overview*

The entire system must satisfy several lofty objectives. The final product must be portable, lightweight, and Internet-enabled. The system must also make measurements of vital bodily functions, perform simple computations, provide data logging functionality, and indicate when measured vitals exceed given ranges, or the user fails to comply with a prescribed logging regimen.

The initial Phase 1 functional requirements for the system are:

- Provide continuous sensor monitoring capability
- Produce visual display of the sensor values
- Accept variety of input data types
- Provide visual indication of warning states
- Provide audible indicator of alarm states

In addition, the following requirements have been added:

- Utilize a hardware-based time reference
- Support dynamic task creation and deletion
- Support a user input device
- Support data logging capabilities
- Support remote communication capability
- Improve overall system performance
- Improve overall system safety

In Phase 3, the following requirements were added:

- Implement a Real time operating system
- Utilize an on-chip thermometer for raw temperature readings
- Display user-selected sensor readings only
- Implement support for EKG sensor readings

In Phase 4, the following requirements were added:

- Implement a web server and web-based user interface

### 3.1.2 Identified Use Cases

Taking the functional requirements listed above, several use cases were developed. A Use case diagram of these scenarios is given in Figure 1. Each use case is expanded and explained below.

#### **Use Case #1: View Vital Measurements**

In the first use case, the user views the basic measurements picked up by the sensors connected to the device.

During normal operation, once the device is turned on by the user, the system records the value output by each sensor. This raw value is linearized and converted into a human-readable form. The user can select toggle between a summary of current vitals as measured by the system or view measurements of each sensor individually.

Three exceptional conditions were identified for this use case:

- *One or more of the expected sensors is not connected* - If this occurs, the measurements taken by the device may be erratic. At the present moment, no action will be taken in such events. Later revisions may address the issue
- *A measured value is outside 5% of the specified normal range* - In this case, a warning signal will flash as an indication of the warning condition
- *A measured value falls outside 10% of a specified "normal" range* - In this case, an audible alarm will sound to indicate the alarm condition

#### **Use Case #2: Acknowledge Alarm**

In the second case, the system is in an alarm state. The user acknowledges the alarm condition by pressing a button.




Figure 1: Use case diagram

Upon pressing the button, the system silences the audible alarm. Any visual warnings continue to flash during the silenced period. If a significant amount of time passes and the sensor reading(s) continue to maintain an alarmed state, the audible alarm will recommence.

Identified exceptions to use case:

- *Alarm is never acknowledged* - If the user never acknowledges an alarm, the system maintains the audible alarm.
- *User continually presses acknowledge button* - The system will only check for a key press at the end of the silencing period. Therefore, the alarm will be silenced if the button is pressed continuously.
- *No speaker or auditory device present* - The auditory alarm is produced via an onboard speaker. If this speaker is removed, no additional audible alarm is supported. This may change in future implementations.

### Use Case #3: View Measurement Logs

In the third use case, the user wishes to view previously recorded measurements for a given vital sign.

As the device is running, the user opts to enter View Logs mode. From this point, the user can select which vital sign they wish to examine. Upon selection, the data logged from previous measurements is displayed.

Possible exception conditions may include the following:

- *User wishes to view more data than the machine remembers* - The machine will not support this operation. The system is only able to display the amount of data defined by the device. Future variants may support an external storage or logging functionality
- *Machine loses power while reading or writing* - if the system loses power, any data stored in dynamic memory will be lost. On restart, data will be overwritten and lost. Future device versions may allow for storage in nonvolatile memory
- *Ongoing measurements trigger warning or alarms* - since measurements are taken continuously, the device may enter an alarm or warning state. In this case, the display will not change unless prompted to by the user. Any alert indicators (visual, audible, or remote messaging) will operate as normal in the background

#### **Use Case #4: System Alert to Remote Terminal**

In the event that the system enters an alert state (e.g. alarm or warning), the system can send a message to a remote terminal connected to the device. The messages sent can inform a second actor of the cause of alert and provide any additional useful information.

Possible exception conditions may include:

- *Improper configuration of the Remote Terminal* - If the remote terminal connection is improperly configured or initialized, data received may be corrupted and not display properly. The device cannot ensure a proper connection and it is the responsibility of the remote user to ensure the correct configuration is used.
- *Remote Connection Lost* - If the user terminates the connection or the connection is lost, messages sent by the device may not arrive at the remote terminal or data may be corrupted. The device will not necessarily monitor the status of any remote connection; this responsibility is the remote terminals. In the event a connection is broken, the device system must continue to perform the other system functions without ill effects.

#### **Use Case # 5: Remote Access via Network**

This device may also be used by an individual at a location outside of normal doctor or medical facilities.

In these cases, a medical professional may be unable to provide direct intervention. Instead, the patient can be monitored remotely from the device through the use of an Internet connection. The remote observer should interact with the device through a terminal just as though they were using the local keypad. Such interactions include starting and stopping measurements, changing measurement selections, and viewing recent measured values.

Possible exception conditions may include:

- *Loss of Network Connection* - The nature of Internet carries the possibility of an unexpected loss of connection. In these instances, information sent from the device or remote terminal may become lost. Various network protocols exist which address and handle these issues. However, the device cannot rely on the remote terminal for operation. In the event data is lost, the local system must continue to operate according to the last received commands or from subsequent local commands.
- *Corruption of Data* - In a remote setting, many environmental variables may affect connection quality. Data sent to and from the device may be corrupted upon reception. The system must be capable of handling corrupted data. In the case of undecipherable commands, the system must fail gracefully.



- *Improperly Formatted Command* - The remote user may send a command that is unsupported by the system. In this case, the user must be informed that an improper command was sent.

### 3.1.3 Detailed Specifications

For this project, the requirements have been further specified as follows:

The system must have the following inputs:

- Alarm acknowledgment capability using a push button
- Buttons or switches to allow user to access system modes and menu items
- Ethernet and webserver to process incoming data requests
- Sensor measurement input capabilities consisting of:
  - \* Body temperature measurement
  - \* Pulse rate measurement signal
  - \* Systolic blood pressure measurement
  - \* Diastolic blood pressure measurement
  - \* EKG frequency measurement

The system must have the following outputs:

- Visual display of the following data in human-readable formats:
  - \* Body temperature
  - \* Pulse rate
  - \* Systolic blood pressure
  - \* Diastolic blood pressure
  - \* Battery status
  - \* EKG Frequency
- Visual indication of warning state with a flashing LED
- Visual indication of a low battery state with an LED
- Audible indication of an alarm state using a speaker
- External data connection to a remote terminal
- Internet connectivity over an Ethernet connection

The initialization values, normal measurement ranges, displayed units, and warning and alarm behaviors for each vital measurement are given in Table 1. The sensors must be sampled every five seconds and the system cannot block and cease operation for five seconds.

A measurement enters a warning state when its value falls outside the stated normal range by 5%.

An alarm state occurs when a measured value falls outside of its specified normal range by more than 15%.

Additionally, the system must be implemented using the Stellaris EKI-LM3S8962 ARM Cortex-M3 microcomputer board, The software for the system must be written in C using the IAR Systems Embedded Workbench/Assembler IDE.

Measurement	Units	Initial Value	Min. Value	Max. Value	Warning Flash Period
Body Temperature	C	75	36.1C	37.8C	1 sec
Systolic BP	mm Hg	80	-	120 mmHg	0.5 sec
Diastolic BP	mm Hg	80	-	80mmHg	0.5 sec
Pulse Rate	BPM	50	60 BPM	100 BPM	2 sec
EKG Frequency	Hz	-	-	-	-
Remaining Battery	%	200	40 %	-	Constant

Table 1: Specifications for measurement data

### 3.1.4 Detailed Task Specifications

- **KeypadTask** The keypad allows the user to interact with the system locally. The user can acknowledge alerts and navigate the menu options.
  - The keypad task will scan the keypad and decode any keypresses
  - The task will have support the following user inputs:
    - Mode selection between 2 modes:
      - \* Measurement Select Menu
      - \* Annunciation
    - Menu selection between 6 options:
      - \* Scan all or measure a specific sensor (Temperature, Blood Pressure, Pulse Rate, or EKG)
    - Alarm acknowledgement
    - Up and down scroll functionality
    - A new set of global variables will be created to store the state of the keypad and key presses
- **Initialize (StartupTask):** The startup task sets up the system hardware when the device is first powered on.
  - The Startup task must be the first task to run
  - It must not be part of the task queue and must only run once
  - The task must configure and activate the system timebase
  - Configure and initialize all hardware subsystems
  - Enable any necessary interrupts
  - Assign priorities to each task prior to creation of the task queue
- **Serial Communication:** Serial communication produces an alert via RS-232 connection to a remote terminal. The connection is strictly uni-directional, serving to alert the remote station that some aspect of the system is in a warning or alarm state.
 

There are no changes to the Serial Communication task since lab 3.

  - The task is enabled by the warn/alarm task

- When run, the task will open an RS-232 connection at 115,200 baud, no flow control, no parity, and 1 stop bit
- The present corrected measurement will be displayed on the terminal in the same fashion as the display task annunciation mode. See Figure 2
- After sending data to the terminal, the serial communication task will remove itself from the task queue




Figure 2: Expected screen layout of remote serial terminal

- **MeasureTask:** The measure function performs measurement collection of the various sensors attached to the system. The current sensors are a thermometer, blood pressure sensor, and heart rate monitor.

No new changes were made to the measure task

- Measurements must be captured every 5 seconds
- Once a complete set of measurements has been taken, the compute task is added to the task queue
- Pointers to the variables used in the measure task will be relocated to accommodate the new data architecture
- The pulse measurement will monitor and count the frequency of a pulse rate event interrupt
- A new value will be stored to memory if the present reading is greater than  $\pm 15\%$  of the previous measurement
- The measurement limits will correspond to 200bpm and 10bpm, determined empirically.

- **ComputeTask:** The compute task performs numerical calculations to linearize and convert the raw sensor measurements into human-readable values.

No new changes were made to the compute task

- The following measurements will be recomputed each time the compute task is run
  1. Temperature

2. Systolic blood pressure
  3. Diastolic blood pressure
  4. Pulse rate
  5. EKG frequency
- This task must only be scheduled when new data is available from the Measure Task
  - The Compute Task must apply any linearization corrections to the raw data
  - The Compute task will store the computed and corrected values into a computed data buffer
  - After computing the corrected values for all measurements, the ComputeTask will remove itself from the task queue
- **DisplayTask:** In the Display task, data collected by the system is displayed on the local OLED display based upon user input received from the keypad. In addition to simply displaying data, the display task will format the data and present a series of menus for the user. The front panel of the OLED display in the Annunciation state is shown in Figure 3a
- No new changes were made to the display task
- Display must support multiple display options
  - Menu mode will allow selection of each of the individual measurements. Upon selection of a measurement, the current value of the measurement will be displayed onscreen
  - Annunciation mode will display the current status of each measurement as in project 1, and provide the same functionality as the display in project 1.
  - The display screen must appear similar to those given in Figure 3
- **Warn/AlarmTask:** The warn/alarm task performs system monitoring functions, creating alerts when measured values exceed the normally accepted ranges. No changes were made to this task in Phases 3 or 4
- The warnings will be activated and indicated as before in project 1
  - The alarm state will be triggered whenever any value is outside 15% of the normal range
  - The alarm will sound in 1 second tones (1 second on, 1 second off)
  - When an alarm or warning state occurs, the serial communication task will be added to the task queue
  - The deactivation period of the alarm sound is defined as 5 measurement periods
- **Schedule:** The Schedule function has been replaced by the RTOS scheduler. Just as before, the Scheduler manages the order in which tasks are executed. It manages resources based on the current running task needs and guarantees tasks complete within the desired timeframe.
- The Scheduler will provide non-preemptive priority-based scheduling

Temperature: 36 C  
Systolic Pressure:  
191 mm Hg  
Diastolic Pressure:  
72 mm Hg  
Pulse rate: 98 BPM  
EKG: 35 Hz  
Battery: 0 %  
IP: 128.95.141.192

(a) Annunciation mode

Make Selection  
->Blood Pressure  
Temperature  
Pulse Rate  
EKG  
Battery  
IP: 128.95.141.192

(b) Menu Select mode

Temperature:  
36 C  
IP: 128.95.141.192

(c) Temperature measurement

Battery:  
92 %  
IP: 128.95.141.192

(d) Battery status

Blood Pressure:  
Systolic: 183 mm Hg  
Diastolic: 75 mm Hg  
IP: 128.95.141.192

(e) Blood pressure measurement

EKG:  
0 Hz  
IP: 128.95.141.192

(f) EKG measurement

Pulse Rate:  
98 BPM  
IP: 128.95.141.192

(g) Pulse rate measurement

Figure 3: Expected OLED display screens

- The hardware timer will provide a system interrupt every 250ms or equal to the minor cycle, whichever is shorter
- Task Control Blocks will interact with other tasks via flags and communication buffers
- Tasks will run at least once every five seconds
- The scheduler cannot block for five seconds
- **StatusTask:** This task monitors the status of the onboard battery.
  - Each time the status task is called, the current amount of battery charge left is queried and stored
- **New Task: EKG Capture** Together EKG Capture and EKG Process will receive EKG sensor information and convert it to human-usable data for analysis
  - The EKG Capture function must convert a time-varying sinusoidal signal to a digital format
  - Each sample must range from 0-3Vdc
  - A total of 256 samples per sequence must be taken at 8-bit precision
  - The capture function must be capable of discerning inputs signals ranging from 35Hz to 3.75kHz
  - Once a sequence is complete, the EKG Process function must be called
- **New Task: EKG Process** The EKG Process function transforms the raw EKg sampled sequence and converts it to a single frequency for analysis.
  - The EKG Process function must perform a Fast Fourier Transform (FFT) on the raw data
  - Following conversion to the frequency domain, the FFT result must be converted to the expected frequency range
  - The 16 most recent EKG results must be stored in a buffer for transmission or review
- **New Task: Command** When interacting with the device remotely, the command task serves as the interpreter of console/browser commands. It conveys system information to the browser in response to webserver requests.
  - Command must be scheduled whenever a command has been received by the web server or when a message must be transmitted to the remote computer.
  - Receive mode:
    - \* A received command must be interpreted and acted upon if valid
    - \* A response acknowledgement or non-acknowledgement must be returned to the webserver whenever a command is received
    - \* Legal commands are given in Appendix B
  - Messages to be sent to a remote computer must be compiled into proper HTML code before being sent to the web server

- **New Task: Remote Communications** Remote communications, or the web server, operates as the gateway between the device and the wider internet. It handles communications protocols and incoming requests to the system.
  - The remote communication task must initialize the network interface
  - This task must also connect to and configure a local area network, handling TCP/IP, HTTP, and other protocols
  - Set up a webserver and webpage to act as a gateway for remote users
  - Update the webpage as needed
  - The user interface must accept a textual command entered by the user as defined in Appendix B

The web browser user interface must contain a text input box and display information in a manner consistent with Figure 4.




Figure 4: Expected Web Browser Interface

## 3.2 Software Implementation

A top-down design approach was used to develop the system. First, a functional decomposition of the problem was carried out based on the identified use cases. Next, the system architecture was developed. After understanding the system architecture, the high-level project file structure in C was defined, followed by the low-level implementation of the tasks.

### 3.2.1 Functional Decomposition

After understanding how the user would interact with the device, the high level functional blocks were developed. These blocks are shown in Figure 5.

The functional blocks were then refined further, showing the main functions the system needed to perform. This refined functional decomposition is shown in Figure 6.






Figure 5: Top level system functionality




Figure 6: Low Level Functional Decomposition

### 3.2.2 System Architecture

Next, the system architecture was developed (Figure 7). At a high level the system works on two main concepts, the scheduler and tasks. In this project, FreeRTOS was used as the scheduler. Tasks embody some sort of work being done, and the scheduler is in charge of determining the speed and order in which the tasks execute. Under FreeRTOS, each task is given a priority, and tasks that are waiting to execute are started based on their priority. The system has several tasks, each with their own specific job. For modularity reasons, each task should have the same public interface and the scheduler should be able to run each task regardless of that specific tasks job or implementation. Thus the task concept is abstracted into a Task Control Block (TCB), and the scheduler maintains a queue of TCBs to run. The TCB abstraction is shown in Figure 7 using inheritance, and the fact that the scheduler has a queue of TCBs is shown with composition. The core functionality of the system was divided into the following eight main tasks:

**Initialization Task** This task Initializes data structures and does system startup-related jobs. This task is not actually scheduled to run, it only executes a single time at system startup.

**Measure Task** The measurement task is in charge of interacting with the blood pressure, temperature, and pulse sensors. Each of these is simulated. The blood pressure and temperature are simulated in the CPU. The task will measure the pulse rate by parsing an externally generated square wave of varying frequency; the pulse rate being proportional to the frequency.

**Compute Task** Compute takes the simulated raw data and converts it to the correct units of measurement. Raw temperature data to Celsius, blood pressure to mm Hg, and pulse rate to BPM.

**EKG Capture** This task is in charge of getting raw EKG data. This data will be simulated. A sample frequency is defined (based on the specified min and max frequencies and the Nyquist



theorem), and an ADC is triggered at the same speed as the sample frequency. The input to the ADC will be a sine wave of variable frequency, from 37 Hz to 3750 Hz. Because of this, the sample frequency chosen was 8000 Hz (more than twice the max frequency).

**EKG Process** This task is in charge of processing raw EKG data. It will run a FFT over the raw EKG data. The results of the FFT operation are then used to determine the maximum frequency in the signal sampled. These data are stored in system memory.

**Keypad Monitor** Keypad will check the keypad for user input. It should provide the user with four keys: two for scrolling, one for selection, and one to go back. The monitor updates the local visual display and updates internal system states in response to mirror any selections or movements by the user

**Display Task** The display task will show a user interface on the Stellaris OLED. The user will interact with the display using a keypad. Under normal operation, a menu will be displayed asking users which measurement they would like to see. If the user presses back while in this menu, they enter annunciate mode which displays all the measurements currently in warning or alarm state.

**Warning/Alarm Task** Under normal operation, this task will light a green LED signifying that everything is OK. If one of the measurements enters a warning state, the task will flash a red LED at a rate specific to the warning for that measurement. If there is an alarm state, it sound the alarm by driving the speaker. At any time if the battery goes too low, the yellow LED will illuminate.

**Serial Task** This task is in charge of sending data to a remote terminal. If any of the states are in a warning or alarm condition, this task will transmit the (corrected) data to the remote terminal. The displayed data will be continually updated until the system returns to a normal condition. Once in a normal condition, the last non-normal state is displayed.

**Command Task** The command task should parse commands from the remote communications task and send a response. The response is based on the particular command sent. For example, if the user requests to view all measurements, the response should be an HTML formatted response containing all of the data for each measurement. Each response should also contain an ACK or NACK, depending on whether or not the command was valid or not.

**Remote Communications** Remote Communications is not actually a task. Instead, the http server, httpd, that is shipped along with the TI embedded development kit was used. This server is interrupt driven and runs entirely in the background. It communicates with tasks using buffers.

**Status Task** Status receives information about the battery on the system and updates its current data accordingly.

Each of these tasks interact using the shared data shown in Figure 7.




Figure 7: System Architecture Diagram

### 3.2.3 High-level Implementation in C

After developing the system architecture, the design needed to be translated into the C programming language. The design manifested in a multi-file program consisting of the following source files:

- **globals.c, globals.h** - Used to define the Shared Data used among the tasks
- **timebase.h** - Defines the timebase used for the scheduler and tasks
- **task.h** - Defines the TCB interface for a task

Each task also has its corresponding “.c” and “.h” file (for example, measure.c and measure.h).

The TCB structure that the scheduler uses must work for all tasks, and must not contain any task-specific information. Instead, the TCB consists of only a void pointer to the tasks data, and a pointer to a function that returns void and takes a void pointer, as shown in Listing 1.

```

1 struct TCB {
2     void *taskDataPtr;
3     void (*taskRunFn)(void *);
4     void *TCB nextTCBPtr;
5     void *TCB prevTCBPtr;
6 }

```

Listing 1: TCB Construct

Leaving out the type information allows the scheduler to pass the task’s data (\*taskDataPtr) into the task’s run function completely unaware of the kind of data the task uses or how the task works.

For increased modularity, the data structure used by each task was not put in the task's header file. Instead, the structure was declared within the task implementation file, and instantiated using a task initialization function. In the header file, a void pointer pointing to the initialized structure is exposed with global scope, as well as the task's run and initialization functions.

### 3.2.4 Task Implementation

The primary task of this project is to implement C code for a medical device on the Stellaris EKI-LM3S8962 and its ARM Coretex A3 processor. The project was started by creating a main file that initializes the variables used in each task and starts the hardware timer then runs into an infinite while loop. Inside the while loop a run method is called. The run method is part of the scheduler. The run method has a runTask flag which determines whether or not anything should actually be run this call. The flag is set to true by the hardware timer's interrupt handler. Once the runTask flag is true the run method will keep track of whether the device is on a minor cycle or a major cycle and run the preform task method of each task. The runTask flag is then set to false so that the tasks will not be executed again until the hardware interrupt has again been triggered. The tasks included in this project are Compute, Measure, Warning, keyPad, OLEDDisplay, Serial, and status. Each task has a public interface of 2 void pointers. One that when initialized by the main method will point to the preform task function, and another that, when initialized, points to a struct containing pointers to the data required by that task. Each task has a task control block(TCB) in the scheduler. This TCB contains pointers to the preform task function and the data for the task and also has fields for pointers to a next TCB and previous TCB. The TCB is used by the scheduler to run the task. The scheduler contains a doubly linked list of TCB objects that uses the TCB next and previous elements to point to the next and previous tasks in the task queue. In this case there are 7 tasks but not always 7 tasks in the list of tasks to run. The compute task is only to be run after the measure task, and the serial task only needs to be run at certain times. When a task is not being used it's TCB is not included in the linked list of tasks. Therefore an updateQueue method was created. This method checks flags set within the tasks that are running and determines if a task that isn't in the list needs to be added, or if a task that is in the list needs to be removed. The scheduler's run task contains a loop that runs through the linked list of TCBs and runs the function pointed to by the TCB with the argument of the data pointer stored in the TCB until the null value pointed to by the last element in the list is reached. After running all tasks the runTask flag is set to false again. The hardware timer will count up until it reaches the number of ms that corresponds to a minor cycle then trigger a hardware interrupt. The interrupt handler sets the runTask flag back to true and the task linked list will be updated, traversed and run again.

Control flow is shown in Figure 8.

Each task has its own unique purpose in the system, and each uses a different part of the global data.

**Measure Task** The Measure task (Figure 9 in Appendix C), deals only with the raw data from the instruments. This task is meant to act in place of the instruments that are unavailable. The task only runs if the scheduler has set the global value is Major Cycle to 1. On a major cycle the measure function either increments the data of each measurement by 1 or 2 or decrements the data by 1 or 2. In the newly improved design a heart rate monitor has been added. The heart rate monitor uses an interrupt handler to count the number of rising edges on an input in a 2




Figure 8: Control Flow Diagram

Major cycle period. The number is then used to calculate the heart rate the sensor is receiving. Additionally the measure task adds the compute task to the task queue after it has run.

**Compute Task** The compute task is very simple. This task will only run after it has been placed into the task queue by the measure task and it will remove itself from the queue after running. It takes the raw data that has been set by the measure task, multiplies by a constant and adds a constant to each piece of raw data to get the corrected data. The compute task then puts the values for the corrected data in memory at the location of the global data pointer. Compute uses every data value except the battery and keypad data. A diagram of the Compute activity is shown in in Figure 10 of Appendix C

**Warning Task** After compute, the warning task begins checking for warning or alarm states. The warning task only deals with the corrected data from compute and the battery state. This task also must deal with the input and output signals used to display warning and sound an alarm. Unlike the other tasks, this task has more to its initialization than just initializing the

data. In addition to setting up the pointers to the global data, during initialization, the task also enables peripheral banks C, F, and G. These are enabled using the `SysCtlPeripheralEnable` library call. Additionally the task set up pins C 5, 6, and 7 as outputs, and pins F 0 and G 1 as PWM outputs. Additionally the PWM outputs are set to use a 65 Hz clock to play a sound at this frequency whenever enabled. The activity diagram is shown in Figure 12 in Appendix C. There are 3 subsystems in the warning task. These subsystems each handle a different part of user notification. The first subsystem deals with the alarm. The subsystem checks to see if the systolic blood pressure data is out of range by more than 20%. If the value is out of that range then the PWM output is enabled using `PWMGenEnable` and pulsed with a 2 second period. If the values fall back within the acceptable range, or the user hits the acknowledge button, then the PWM is disabled with `PWMGenDisable` and the sound stops. The acknowledge button is sensed in the keyPad task and a global data value `global.alarmAcknowledge` is set. When the `alarmAcknowledge` field is set to 1 the alarm will go to its resting state. The next subsystem checks the corrected data for being 5% out of range of the accepted values. If any value is more than 5% out of its range then a warning will be displayed on the red led connected to pin C 5 using the `GPIOPinWrite` function. Depending on the value that is out of range the period the led flashes at will vary. In addition to the led flashing the warning will also tell the scheduler to add the `serialCommunication` task to the task list. The final subsystem is the battery check. This system checks if there is more than 30% of battery left on the device. This is taken from the battery state data field. If there is less than 20% battery left then a yellow led connected to pin C7 is illuminated, if there is more then 20% battery, and the device is not in a warning or alarm state then the green led on pin C 6 is illuminated.

**Keypad Task** The keyPad task uses the GPIO libraries to set up 5 inputs. 1 input is pin E 0 and is used as the alarm acknowledge button. The other 4 pins are inputs on pins D 4, D 5, D 6, and D 7, and are used to detect button presses on the external keypad. The keypad works by connecting 2 of the outputs from the keypad together. There are 4 row outputs and 4 column outputs for a total of 16 possible combinations. The design was simplified by only using 1 column of 4 buttons. By doing this, the column keyPad wire can be supplied with a constant 5V from the Stellaris board, then each of the 4 row lines can be connected to a GPIO input. when a button in the live column is pushed the 5V column line is connected to 1 of the GPIO pins which lets the keyPad task know which button was pushed. Every time the keyPad task is run the GPIO will tell the task what buttons are pressed. 2 of the buttons correspond to up and down for scroll in the menu, 1 button corresponds to select in the menu, and 1 button changes the mode between menu and annunciation. The data that the buttons manipulate are the global keyPad data parts which are `global.select`, `global.scroll`, and `global.mode`. Additionally the input on pin E 0 manipulates the `global.alarmAcknowledge` signal. The activity diagram for this task is shown in Figure 11 in Appendix C.

**Display Task** To show a user their current medical measurements, the system also has an `oledDisplay` task. This task uses the corrected data from measurements, the battery state, and the keyPad data. The display task has an activity diagram shown in Figure 14 in Appendix C. This task uses the `usnprintf()` function in C to convert the data types that the corrected data is stored in, and properly format these data values into a string which is stored in a buffer. The string contained in the buffer is then printed to the OLED screen using the driver library `rit128x96x4` functions. The display has 2 modes. The mode that the screen currently displays is determined by the `global.mode` data which is set in the keyPad task. When mode is 0 this is the menu mode. This mode displays each of the 4 measurement types, temperature, blood

pressure, pulse rate, and battery without the data for each. The task then takes the global.scroll data from keyPad and displays a cursor next to the measurement that scroll is currently at, 0 corresponding to temperature, 1 to blood pressure, and so on. If the global.select is set to 1 in keyPad then the currently scrolled to measurement is selected and the screen will change to showing only that measurement and its data. When global.mode is equal to 0 the annunciation screen will instead be displayed. This screen shows each measurement followed by its data.

**Serial Task** The serial task is only run when a warning occurs. The warning causes the TCB for the serial task to be added to the TCB schedule linked list. When the serial task is run the first time it will initialize a UART connection using UARTConfigSet and UARTEnable driver functions to enable the UART that communicates to an FTDI chip on the Stellaris board. The FTDI chip then converts the UART to a virtual serial port over the USB cable to the PC. After the UART is initialized on the first run and in all subsequent runs of the serial task, all the measurements and their data are formatted and printed into 1 buffer using the usnprintf function. A loop then iterates through the buffer writing each character in the buffer to the UART. An activity diagram showing the serial task is located in Appendix C as Figure 13.

**Status Task** The last task is the status task. This task only deals with one piece of data which is the battery state data. The only thing the status task does is that on a major cycle, it decrements the battery state by 1. This is shown in the activity diagram in Figure 15 in Appendix C.

## 4 PRESENTATION, DISCUSSION, AND ANALYSIS OF THE RESULTS

### 4.1 Results

The project was completed and demonstrated on March 17, 2014.

Demonstration of the system to the interested parties showed that the system met the majority of the requirements initially presented at the onset of the lab project. Testing of the system prior to demonstration also verified that the system met the specifications listed in Section 3.1.

During the demonstration, all tasks worked as designed and expected with the exception of the serial communication task, serialTask. After initially performing as expected, the system would freeze up. When this happened, the system was unresponsive and did not produce any annunciation.

Using an oscilloscope, the run times of each task were empirically determined. The procedure used involved asserting a GPIO pin high immediately after entering the function in question and deasserting the pin just prior to exiting the function. The trace of the pin output helped determine the exact function timing. The results are given in Table 2

**Answers to the last three questions in the list of items to include in the project report:**

*You don't find the stealth submarine. That's why they are so expensive; at that cost, you take great pains to never lose one.*

*A helium balloon always rises. It just rises upside-down.*

*If you really managed to lose the stealth submersible, you first have to tell the government, which will deny it has any stealth submersibles, then you have to comb the seven seas until your comb hits the sub.*

Task	Runtime ( $\mu$ s)
Measure	20.3
Compute	55.4
Display	33200
Warning	17.7
Serial	12000
keypad	5.72
ekgCapture	26600
ekgProcess	1345
Status	5.6

Table 2: Empirically determined task runtimes

## 4.2 Discussion of Results

The ease of change in the code is the result of a large amount of time spent on design. The design makes it easy to configure flash times, add new tasks, and to reason about tasks independently of the whole system. The solid high-level architectural design led to ease of implementation and change.

In terms of performance, the run times of each task appear to correspond with the number of instructions required for each task. Given the speed of the CPU, 8 MHz, we can calculate an estimated number of instructions for each task. This is given in Table 3. The majority

Task	Instructions
Measure	162
Compute	443
Display	265600
Warning	142
Serial	96,000
keypad	45.8
ekgCapture	212,800
ekgProcess	10,760
Status	45

Table 3: Estimated instructions per task, rounded to the nearest instruction

of the cycles are likely spent waiting for memory. For example, the status task only has two comparisons and an arithmetic operation, but has to reference the data in global memory. The exception here is the display task, which was about three orders of magnitude more instructions than the other tasks. This was due to the `sprintf()` library call, included in the standard C library. While this could have been optimized, it was found that with a minor cycle delay of 250 ms, the display delay of 33.2 ms was not significant.

## 4.3 Analysis of Any Errors

There were two errors found in the final project. Both errors involved a lack of accuracy in the reported values measured by the sensors. The pulse rate measurement was off in the same manner as in Lab 3. The other measurement, from the EKG sensor, was off by a similarly small amount. The cause of this error differed however.

A previously discussed error in the serial terminal display has been successfully resolved.

The specified corrected pulse rate was to be between 10 BPM and 200 BPM. For simplicity, the raw pulse rate was implemented as a 1-1 mapping from frequency to raw pulse rate. For example, a 1 Hz frequency would produce a value of 1 for raw pulse rate, and a 15 Hz signal would produce a value of 15 for raw pulse rate. This caused two issues. When the Input frequency is 1, the measured BPM (using the specified raw to correct conversion) is  $8 + 3 \cdot (1) = 11BPM$ . This is larger than the specified 10 BPM minimum. Also, when the input frequency is 64 Hz, a corrected value of  $8 + 3 \cdot (64) = 200BPM$  is expected. However, as the frequency increased, the overhead of the other running tasks became significant. As a result, more rising edges could fit in our measurement interval than we expected. This resulted in a maximum BPM of roughly 206 BPM.

In the case of the EKG measurements, the accuracy of the reported value may vary by as much as  $\pm 5\%$  BPM. The actual inaccuracy is a function of the actual EKG frequency. We have attributed this inaccuracy to the implementation of the EKG measurement functions. Specifically the timing mechanism used sample the analog sensor relies on a software defined delay function to capture the signal. This choice was made because the onboard hardware did not have a sufficient number of hardware based timers to use. As a result, the EKG measurements are forced to rely on a system that is not as accurate initially. Previous experience has indicated that a software timer can made highly accurate given enough calibration and adjustment of correcting variables. Time constraints limited our calibration to a lower level of accuracy for the time being. In addition, a more precise determination of the actual sampling rate will yield correspondingly more accurate measurements of the EKG frequency. We have made a decision to consider the level of precision sufficiently accurate for this prototype version. A subsequent version will build upon these lessons learned here.

#### **4.4 Analysis of Implementation Issues and Workarounds**

The medical instrument design in this project was completed and tested successfully to meet almost all the requirements, the designers did face a few difficulties in designing the device, however, because this design was additional functionality added to a previous projects design, many of the errors previously encountered were easily avoided due to experience of the students, and already completed coding work.

Many of the challenges the designers of this project faced were in the keypad input and the data output. The keypad input posed a difficult hardware challenge as there are 16 input keys on the keypad but only 8 connections for the microprocessor to connect to the keypad. This means that to identify a single key press 4 connections must be set as outputs and 4 as inputs. The inputs can then be scanned as the outputs are set 1 at a time to find which key is pressed. Instead of implementing this design, the students instead opted to use only 4 buttons on the keypad. This allowed the strobe design to be ignored. Instead 1 row of keys was activated all the time and that row was scanned for button presses.

In addition to keypad input, there was also difficulty in implementing data formatting functions. After adding a hardware delay, IAR workbench no longer allows the use of `sprintf` which had previously been used to print and format data. The `usnprintf` command was instead used to format and print data to a buffer, however, the students found that `usnprintf` does not have the ability to print floating point data. This issue was resolved by changing measurements that were previously printed as floats to be printed as integers. `usnprintf` also caused issues when printing certain data for the serial task. In this case the `usnprintf` was causing a runtime error and freezing the operation of our device. This issue remained unresolved.



All problems but one were solved before demonstrating the product to the interested parties. The final project still contained the serial error previously mentioned.

## 5 TEST PLAN

To ensure that this project meets the specifications listed in section 3.1, the following parts of the system must be tested:

### **Phase I Tests:**

- Vitals are measured and updated
- System properly displays corrected measurements and units properly
- System enters, indicates, and exits the proper warning state for blood pressure, temperature, pulse, and battery
- System enters and exits the alarm state correctly
- Alarm is silenced upon button push
- Alarm recommences sound after silencing if system remains in alarm state longer than silence period

Additional tests to determine the runtime of each specific task are also required.

The inclusion of additional specifications for Phase II of the project requires additional tests to ensure the system meets the customer requirements.

### **Phase II Tests:**

- Scheduler loops through linked list properly
- Scheduler adds and removes from the linked list the following tasks correctly:
  - Compute task added by measure task
  - Serial communication task added by annunciation task
  - Compute task removed by itself
  - Serial communication task removed by itself
- Warning task alarm meets the following two requirements
  1. Has one (1) second tones; a total period of 2 seconds
  2. Activates only when systolic pressure is greater than 20% above normal
  3. Has an auditory deactivation or “sleep” period of 5 measurement cycles
- Serial task displays the temperature, blood pressure, pulse rate, and battery status as listed in Section 3.1
- Keypad task captures user inputs, sets the appropriate inputs, and causes the associated changes in system state
- Hardware timer updates the system’s minor cycle counter

The specifications added during Phase III require that the following Phase III properties be tested to satisfaction.

**Phase III Tests:**

- EKG related functions are able to accurately and reliably measure signals within the specified ranges
- The System operates using an non-preemptive, priority-based real time operating system
- Users may use the keypad to retrieve and view specific measurements
- The system startup task properly initializes the hardware needed for operation

The following tests are required to ensure the specifications and requirements for Phase IV are met.

**Phase IV Tests:**

- Command task properly handles expected commands
- Command Task properly handles unexpected and incorrect commands
- Web server properly responds to browser requests and forwards information to the system
- Remote connection can initialize and maintains a connection for a significant period of time
- User interaction from keypad and web browser can occur together without significant loss of functionality

More detailed explanation of the tests performed is provided in the following sections.

## **5.1 Test Specification**

### *5.1.1 Scheduler*

The scheduler (FreeRTOS) needs to be shown to correctly schedule and dispatch tasks. This means that task should execute in the right order, and at the right time. Given a minor cycle of 50 ms, every task should run roughly once every 50 ms. Also, the scheduler needs to successfully add and remove tasks from the queue dynamically. Specifically, the Measure Task should be able to tell the scheduler to add the Compute Task and the Warning/Alarm Task should be able to schedule the Serial Task. Both Compute and Serial should be able to be removed from the schedule.

### *5.1.2 Measure Task*

For this design, the temperature and blood pressure values were simulated on the CPU. The pulse rate was simulated using an externally generated square wave of varying frequency.

- **Temperature** The temperature should increase by two every even major cycle (5 seconds) and decrease by one every odd major cycle until it exceeds 50, at which point the process should reverse (decrease by two every even major cycle and increase by one every odd major cycle), until it dips below 15, and the whole process should be started over again.

- **Pulse** The pulse rate should match one-to-one with the frequency of the input signal. For example, a 15 Hz signal should produce a raw pulse rate of 15.
- **Systolic Pressure** The systolic pressure should increase by three every even major cycle and decreases by one every odd major cycle. If it exceeds 100, it should reset to an initial value.
- **Diastolic Pressure** The diastolic pressure should decrease by two on even major cycles and decrease by one on odd major cycles, until it drops below 40, when it should restart the process.

The Measurement Task should also successfully add the Compute Task to the schedule queue.

### 5.1.3 Compute Task

The compute task should be verified to convert raw simulated sensor data according to the following formulas.

- $CorrectedTemperature = 5 + 0.75 * RawTemperature$
- $CorrectedSystolicPressure = 9 + 2 * RawSystolicPressure$
- $CorrectedDiastolicPressure = 6 + 1.5 * RawTemperature$
- $CorrectedPulseRate = 8 + 3 * RawTemperature$

The compute task should also successfully remove itself from the schedule queue.

### 5.1.4 Keypad Task

The keypad should be tested to successfully capture user input. When the select button is pressed, the measurement selection value should reflect the selected measurement. When the up scroll button is pressed, the scroll value should be incremented, and when the down scroll button is pressed, the scroll value should be decremented. If the alarm acknowledge button is pressed, this should be reflected in the alarmAcknowledge global value.

### 5.1.5 Display Task

On load, the display task should present the user with an option to select the desired measurement. If the back button is pressed, the annunciation screen should be displayed, showing the measurements in warning or alarm state.

### 5.1.6 Warning/Alarm Task

The warning/alarm system needs to be tested to do several things. When in a warning state, it should flash the red LED at the rate appropriate for the warning. When the battery is low, it should illuminate the yellow LED. If the system is in an alarm state, it should sound the speaker alarm. The following ranges in Table 4 are calculated from the specified minimum and maximums found in Table 1 on page 6.

This task should also add the Serial task if any of the measurements are in a warning or alarm condition.

Data	Warning Range	Alarm Range
Temperature	34.3 - 39.7 C	32.5 - 41.6 C
Systolic Pressure	> 84 mmHg	> 88 mmHg
Diastolic Pressure	> 126 mmHg	> 132 mmHg
Pulse	57 - 63 BPM	54 - 110 BPM

Table 4: Initial values and warning/alarm states

#### 5.1.7 Serial Task

If any of the measurements are in warning or alarm state, this task should send this data serially to a remote terminal. The task should send all the data (not just the data in warning or alarm state). It should be printed as shown in Listing 2.

```

1 1. Temperature          0 C
2 2. Systolic Blood Pressure 0 mm Hg
3 3. Diastolic Blood Pressure 0 mm Hg
4 4. Pulse Rate          0 BPM
5 5. Battery              0 %

```

Listing 2: Remote Terminal Output

#### 5.1.8 Status Task

Since the initial design does not use a battery, the status task simulates the battery state using the CPU. For now, it simply decrements the state of the battery. The test should show that the battery state is decremented by one every major cycle.

#### 5.1.9 EKG Capture

EKG Capture should be tested to grab correct values from an input source (sine wave).

#### 5.1.10 EKG Process

EKG Process should successfully measure sine waves between 37 and 3750 Hz. Key frequencies should be chosen between this range, like 37, 40, 60, 100, 200, 500, 1000, 3000, and 3600, and 3750. The system should successfully detect each of these frequencies.

#### 5.1.11 Command Task

The command task should successfully parse each of the commands shown in Appendix B. Upon receiving a command, it should also generate a correct response. For example, a “M T” command should give the response shown in Listing 3. The “A” refers to an ACK. If the command was invalid, the response should simply be an “E”.

```

1 A
2 <p>Temperature 0 C</p>

```

Listing 3: Command Response

### 5.1.12 Remote Communications

The server should be tested to successfully respond to HTTP requests. This means opening a file and sending it over HTTP to the requester. Since the system also has special URLs for commands, it should also recognize these commands.

## 5.2 Test Cases

The students begin testing by examining if the alarm sounds at the proper time. This is initially tested by disabling the functions that simulate measurements being made on each of the data measurements, and setting their initial values to be either within the alarm range or outside of the alarm range. The warning states were also initially tested this way. The initial values for raw data given in Table 1 on page 6 were used to test the normal state of the machine because each falls within the acceptable range of measurements for corrected data (also given in Table 1) that does not require a warning. Using these initial values, the code was programmed onto the Stellaris board. Correct operation was verified by the alarm not sounding, and the red led being off, indicating that no warning state was in effect. In addition the green led was on indicating a normal state. Next the students varied one parameter at a time to be outside of the acceptable range by more than 10%. Starting with the temperature being set to an initial raw value of 50, the alarm was verified by hearing the aural annunciation coming from the system. In addition, the temperature warning stat was also in effect. This means that the green led was off and the red led was blinking. To verify correct operation we needed to make sure the led was blinking with a period of 1 second. The correct flashing pattern was verified by counting the number of times the led flashed in 6 seconds. In this case, for temperature, the led flashed 6 times in 6 seconds indicating a 1 second period, and correct operation. After this test, the temperature value was returned to 42 and the Pulse was instead set to 45. The same methods were used to verify that the alarm and warning states for pulse rate were working correctly, but this time the warning led turned on 3 times in 6 seconds indicating a 2 second period which is the intended period of flashing. The pulse rate was then returned to 25 and each pressure reading was checked for correct operation individually by being set to an initial raw value of 100. Once again, the green led started off because the system was not in a normal state. The alarm was sounding due to the extremely high blood pressure measurements, and the red warning led flashed 12 times in 6 seconds indicating the correct period of .5 seconds for a blood pressure warning. In addition to testing the validity of each warning state and alarm state, the acknowledgment of the alarm was also tested during each of these tests. This was tested by hitting the acknowledge button once during each measurements test. During each test, hitting the acknowledge button turned the alarm sound off for a short time, as intended.

Next the measurement simulation functions were tested. This was done by re-enabling each one that had been disabled from the previous test one at a time. The initial raw values were again set to the values in Figure 5. When each measurement was re-enabled, the students could watch the temperature change at each major cycle using the OLED display. Since the OLED display indicated that the corrected temperature went up .75 degrees on a major cycle then down 1.5 degrees on the next, the temperature measurement was working as intended. This situation also gave the students an opportunity to verify that the warning and alarm states initiated as the temperature fell out of the acceptable range. The Led began flashing with a 1 second period after a few major cycles, then the alarm began sounding, indicating correct operation. Since temperature was working correctly, the temperature measurement function was once again disabled and each blood pressure measurement was re-enabled individually for testing. The Systolic pressure began by rising 4 mm Hg on a major cycle then falling 2 mm

Hg on the next, and the Diastolic pressure by rising 3 on a Major cycle and falling 1.5 on the next, this was consistent with the design. The warning and alarm states were activated as each passed its threshold and the red led was blinking with a period of .5 seconds. The warning led was also tested in the case that all warning states were active. To do this all initial values were set to 100. In this case, as designed by the students, the red warning led indicated the fastest blinking warning with a .5 second period.

In this device, a new pulse rate monitor device was added. To test the pulse rate monitor the monitor input was connected to a function generator generating square waves. As the square wave frequency was increased the pulse rate value was expected to increase. This was verified to be working correctly. As the pulse rate passed through the warning zones the design for pulse rate warnings was also verified to be working correctly.

Additionally, the improved medical device now has a menu that is displayed on the OLED display and navigated using the testbench keypad. The design operating these functions was tested by navigating to each part of the menu using the keys on the keypad and visually verifying that each part of the menu displayed the correct data. Each menu, annunciation, main menu, and each measurements selection mode, was visually verified to be working as intended. The keypad functionality was verified as each button was used to navigate through the menu.

The newly added serial communication was then tested using the hyperterm program on the lab station PC. The serial connection was established over a virtual COM port on the USB connection from the PC to the Stellaris board. The program had the intended functionality of displaying each of the measurements, and its current data on the serial port whenever the alarm state was entered. The functionality could be verified by watching the hyperterm screen to see if the data displayed when a warning state occurred. The data on hyperterm could then be compared to the OLED display data to verify its accuracy.

The final bit of testing performed on the system was timing each task within the system. This was done by adding a general purpose output pin in the scheduler code. This output was set high right before the execution of a task, and set low immediately after the execution of the task. An oscilloscope was then attached to this output pin and set to trigger on a positive edge. The cursors were then used to measure the amount of time the signal was high in each cycle.

## **6 SUMMARY AND CONCLUSION**

### **6.1 Final Summary**

A medical monitoring system with a user interface, alarm notification, and remote terminal display was designed, implemented, and tested. The design simulated temperature and blood pressure, and obtained pulse rate data from an external function generator. These results were converted to a human readable form, and tested to see if there was a warning or alarm state. In the event of warning or alarm, the user was notified visually with LEDs and aurally with an alarm sound. The warning and alarm data was transferred to a remote terminal.

Some implementation errors were encountered. The pulse rate range could not fit the specification exactly, and the serial communications task was not printing the correct values to the remote terminal. Aside from these two implementation errors, the device worked as specified.

### **6.2 Project Conclusions**

This project contained 3 major phases, the design, implementation, and testing steps. The students were immediately introduced to using the unified modeling language(UML) to design embedded systems. This is the first time many students will have used UML for system design

which caused some confusion and difficulty. In the end through the use of the UML guidelines for design, the students were able to implement their system in code for the Texas Instruments Stellaris EKI-LM3S8962 much more quickly and with far fewer errors than if they had spent less time in the design phase of this project.

Effective design tools allowed the students to quickly implement their embedded system in C code for an ARM Cortex A3 processor, and move onto the testing phase of the project quickly. Unfortunately, while testing the students encountered a number of problems in using the PWM and general purpose input and output signals. After consulting the documentation for the Stellaris kit and solving their input/output problems, they began testing their design using visual and audio cues, the IAR embedded workbench debugger, and a few specifically programmed debug features. After the results of the testing verified the design to be working correctly, the students proceeded to present their medical instrument to their instructor.

**A BREAKDOWN OF LAB PERSON-HOURS (ESTIMATED)**

Person	Design Hrs	Code Hrs	Test/Debug Hrs	Documentation Hrs
Patrick	20	20	30	12+
Jarrett	9	30	16	9
Jonathan	15	30	15	10

By initializing/signing above, I attest that I did in fact work the estimated number of hours stated. I also attest, under penalty of shame, that the work produced during the lab and contained herein is actually my own (as far as I know to be true). If special considerations or dispensations are due others or myself, I have indicated them below.



## **B SUPPORTED WEB BROWSER COMMANDS**

The Commands and Responses for the embedded application task are given as follows.

**S** The S command indicates START mode. The command shall start the embedded tasks by directing the hardware to initiate all the measurement tasks. In doing so, the command shall enable all the interrupts.

**P** The P command indicates STOP mode. This command shall stop the embedded tasks by terminating any running measurement tasks. Such an action shall disable any data collecting interrupts.

**D** The D command enables or disables the OLED display.

**M;payload;** The M command. The M command requests the return of the most recent value(s) of the specified data. The M response. The M response returns the most recent value(s) of the specified data.

The ;payload; is defined as one of the following:

**A** - Instruct the system to measure all sensors. Will return the most recent measurement of each sensor

**T** - Instruct system to measure Temperature. Will return the temperature in degrees Celsius

**B** - Instruct system to measure the patient blood pressure. Will return the systolic and diastolic blood pressure in mmHg

**P** - Instruct system to measure patient pulse rate. Will return the heart rate in beats per minute

**E** - Instruct the system to perform an EKG measurement. Will return the primary EKG frequency

**D** - DEBUG mode only. Returns the local display status

**M** - DEBUG mode only. Returns system measurement status. Returns boolean true if measurement is enabled, false otherwise

**A** The A response acknowledges the receipt of the identified command.

**E** The E error response is given for incorrect commands or non-existent commands.

## C ACTIVITY DIAGRAMS




Figure 9: Measure Activity Diagram




Figure 10: Compute Activity Diagram




Figure 11: Keypad Activity Diagram




Figure 12: Warning Activity Diagram




Figure 13: Serial Activity Diagram




Figure 14: Display Activity Diagram




Figure 15: Status Activity Diagram




Figure 16: EKG Capture Activity Diagram




Figure 17: EKG Process Activity Diagram




Figure 18: Command Activity Diagram




Figure 19: Remote Activity Diagram



## D SOURCE CODE

Source code for this project is provided below.

### D.1 Main Function

../code/main.c

```
1  /*
2     FreeRTOS V7.0.1 – Copyright (C) 2011 Real Time Engineers Ltd.
3
4
5     *****
6
7     *
8     *   FreeRTOS tutorial books are available in pdf and paperback.
9     *   Complete, revised, and edited pdf reference manuals are also
10    *   available.
11    *
12    *   Purchasing FreeRTOS documentation will not only help you, by
13    *   ensuring you get running as quickly as possible and with an
14    *   in-depth knowledge of how to use FreeRTOS, it will also help
15    *   the FreeRTOS project to continue with its mission of providing
16    *   professional grade, cross platform, de facto standard solutions
17    *   for microcontrollers – completely free of charge!
18    *
19    *   >>> See http://www.FreeRTOS.org/Documentation for details. <<<
20    *
21    *   Thank you for using FreeRTOS, and thank you for your support!
22    *
23    *****
24
25    This file is part of the FreeRTOS distribution and has been modified to
```

```

26 demonstrate three simple tasks running.
27
28 FreeRTOS is free software; you can redistribute it and/or modify it
29 under
30 the terms of the GNU General Public License (version 2) as published by
31 the
32 Free Software Foundation AND MODIFIED BY the FreeRTOS exception.
33 >>>NOTE<<< The modification to the GPL is included to allow you to
34 distribute a combined work that includes FreeRTOS without being obliged
35 to
36 provide the source code for proprietary components outside of the
37 FreeRTOS
38 kernel. FreeRTOS is distributed in the hope that it will be useful, but
39 WITHOUT ANY WARRANTY; without even the implied warranty of
40 MERCHANTABILITY
41 or FITNESS FOR A PARTICULAR PURPOSE. See the GNU General Public License
42 for
43 more details. You should have received a copy of the GNU General Public
44 License and the FreeRTOS license exception along with FreeRTOS; if not
45 it
46 can be viewed here: http://www.freertos.org/a00114.html and also
47 obtained
48 by writing to Richard Barry, contact details for whom are available on
49 the
50 FreeRTOS WEB site.
51
52 1 tab == 4 spaces!
53
54 http://www.FreeRTOS.org – Documentation, latest information, license and
55 contact details.
56
57 http://www.SafeRTOS.com – A version that is certified for use in safety
58 critical systems.
59
60 http://www.OpenRTOS.com – Commercial support, development, porting,
61 licensing and training services.
62 */
63
64 /*
65 * Creates all the application tasks, then starts the scheduler. The WEB
66 * documentation provides more details of the standard demo application
67 * tasks.
68 * In addition to the standard demo tasks, the following tasks and tests are
69 * defined and/or created within this file:
70 *
71 * "OLED" task – the OLED task is a 'gatekeeper' task. It is the only task
72 * that
73 * is permitted to access the display directly. Other tasks wishing to
74 * write a
75 * message to the OLED send the message on a queue to the OLED task instead
76 * of

```

```

65  * accessing the OLED themselves. The OLED task just blocks on the queue
    waiting
66  * for messages – waking and displaying the messages as they arrive.
67  *
68  * "Check" hook – This only executes every five seconds from the tick hook.
69  * Its main function is to check that all the standard demo tasks are still
70  * operational. Should any unexpected behaviour within a demo task be
    discovered
71  * the tick hook will write an error to the OLED (via the OLED task). If
    all the
72  * demo tasks are executing with their expected behaviour then the check
    task
73  * writes PASS to the OLED (again via the OLED task), as described above.
74  *
75  * "uIP" task – This is the task that handles the uIP stack. All TCP/IP
76  * processing is performed in this task.
77  */
78
79
80
81
82  /*****
83  * Please ensure to read http://www.freertos.org/portlm3sx965.html
84  * which provides information on configuring and running this demo for the
85  * various Luminary Micro EKs.
86  *****/
87
88  /* Set the following option to 1 to include the WEB server in the build. By
89  default the WEB server is excluded to keep the compiled code size under the
    32K
90  limit imposed by the KickStart version of the IAR compiler. The graphics
91  libraries take up a lot of ROM space, hence including the graphics libraries
92  and the TCP/IP stack together cannot be accommodated with the 32K size limit
    . */
93
94  // set this value to non 0 to include the web server
95
96  #define mainINCLUDE_WEB_SERVER    0
97
98
99  /* Standard includes. */
100 #include <stdio.h>
101
102 /* Scheduler includes. */
103 #include "FreeRTOS.h"
104 #include "task.h"
105 #include "queue.h"
106 #include "semphr.h"
107
108 /* Hardware library includes. */
109 #include "hw_memmap.h"
110 #include "hw_types.h"
111 #include "hw_sysctl.h"

```

```

112 #include "sysctl.h"
113 #include "gpio.h"
114 // #include "glib.h"
115 #include "rit128x96x4.h"
116
117 /* Demo app includes. */
118 // #include "lcd_message.h"
119 // #include "bitmap.h"
120
121
122 #include "timebase.h"
123 #include "schedule.h"
124
125 /* Tasks */
126 #include "tcb.h"
127 #include "startup.h"
128 #include "measure.h"
129 #include "compute.h"
130 #include "display.h"
131 #include "warning.h"
132 #include "keypad.h"
133 #include "serial.h"
134 #include "status.h"
135 #include "keyPad.h"
136 #include "ekgCapture.h"
137 #include "ekgProcess.h"
138 #include "commandTask.h"
139
140 /*-----*/
141
142 /*
143  The time between cycles of the 'check' functionality (defined within the
144  tick hook.
145  */
146 #define mainCHECK_DELAY ( ( portTickType ) 5000 / portTICK_RATE_MS )
147
148 // Size of the stack allocated to the uIP task.
149 #define mainBASIC_WEB_STACK_SIZE ( configMINIMAL_STACK_SIZE * 3 )
150
151 // The OLED task uses the sprintf function so requires a little more stack
152 // too.
153 #define mainOLED_TASK_STACK_SIZE ( configMINIMAL_STACK_SIZE + 50 )
154
155 // Task priorities.
156 #define mainQUEUE_POLL_PRIORITY ( tskIDLE_PRIORITY + 2 )
157 #define mainCHECK_TASK_PRIORITY ( tskIDLE_PRIORITY + 3 )
158 #define mainSEM_TEST_PRIORITY ( tskIDLE_PRIORITY + 1 )
159 #define mainBLOCK_Q_PRIORITY ( tskIDLE_PRIORITY + 2 )
160 #define mainCREATOR_TASK_PRIORITY ( tskIDLE_PRIORITY + 3 )
161 #define mainINTEGER_TASK_PRIORITY ( tskIDLE_PRIORITY )
162 #define mainGEN_QUEUE_TASK_PRIORITY ( tskIDLE_PRIORITY )
163

```

```

164 // The maximum number of messages that can be waiting for display at any
    one time.
165 #define mainOLED_QUEUE_SIZE          ( 3 )
166
167 // Dimensions the buffer into which the jitter time is written.
168 #define mainMAX_MSG_LEN              25
169
170 /*
171  The period of the system clock in nano seconds. This is used to calculate
172  the jitter time in nano seconds.
173 */
174
175 #define mainNS_PER_CLOCK ( ( unsigned portLONG ) ( ( 1.0 / ( double )
    configCPU_CLOCK_HZ ) * 1000000000.0 ) )
176
177
178 // Constants used when writing strings to the display.
179
180 #define mainCHARACTER_HEIGHT          ( 9 )
181 #define mainMAX_ROWS_128              ( mainCHARACTER_HEIGHT * 14 )
182 #define mainMAX_ROWS_96               ( mainCHARACTER_HEIGHT * 10 )
183 #define mainMAX_ROWS_64               ( mainCHARACTER_HEIGHT * 7 )
184 #define mainFULL_SCALE                 ( 15 )
185 #define ulSSI_FREQUENCY                ( 3500000UL )
186
187 /*-----*/
188
189 /*
190  * Configure the hardware .
191  */
192 static void prvSetupHardware( void );
193
194 /*
195  * Hook functions that can get called by the kernel.
196  */
197 void vApplicationStackOverflowHook( xTaskHandle *pxTask, signed portCHAR *
    pcTaskName );
198 void vApplicationTickHook( void );
199
200 /* The tasks */
201 void measure( void *vParameters );
202 void compute( void *vParameters );
203 void display( void *vParameters );
204 void keyPad( void *vParameters );
205 void warning( void *vParameters );
206 void serial( void *vParameters );
207 void status( void *vParameters );
208 void ekgCapture( void *vParameters );
209 void ekgProcess( void *vParameters );
210
211 void command( void *vParameters );
212
213 /*-----*/

```

```

214 |
215 | /*
216 |   The queue used to send messages to the OLED task.
217 | */
218 | xQueueHandle xOLEDQueue;
219 |
220 | /*-----*/
221 |
222 | xTaskHandle computeHandle;
223 | xTaskHandle serialHandle;
224 | xTaskHandle measureHandle;
225 | xTaskHandle displayHandle;
226 | xTaskHandle ekgCaptureHandle;
227 | xTaskHandle ekgProcessHandle;
228 | xTaskHandle commandHandle;
229 |
230 | /*****
231 |  * Please ensure to read http://www.freertos.org/portlm3sx965.html
232 |  * which provides information on configuring and running this demo for the
233 |  * various Luminary Micro EKs.
234 |  *****/
235 |
236 | int main( void )
237 | {
238 |     prvSetupHardware();
239 |
240 |     /* Startup task */
241 |     startup();
242 |
243 |     /* Start the tasks */
244 |     xTaskCreate(measure, "measure task", 100, NULL, 2, &measureHandle);
245 |     xTaskCreate(compute, "compute task", 100, NULL, 3, &computeHandle);
246 |     xTaskCreate(ekgCapture, "ekgCapture task", 1024, NULL, 1, &
ekgCaptureHandle);
247 |     xTaskCreate(ekgProcess, "ekgProcess task", 1024, NULL, 3, &
ekgProcessHandle);
248 |     xTaskCreate(display, "display task", 100, NULL, 5, &displayHandle);
249 |     xTaskCreate(keyPad, "keyPad task", 100, NULL, 4, NULL);
250 |     xTaskCreate(warning, "warning task", 100, NULL, 5, NULL);
251 |     xTaskCreate(serial, "serial task", 100, NULL, 5, &serialHandle);
252 |     xTaskCreate(status, "status task", 100, NULL, 1, NULL);
253 |     xTaskCreate(command, "command task", 100, NULL, 3, &commandHandle);
254 |
255 |     vTaskSuspend(measureHandle);
256 |     vTaskSuspend(computeHandle);
257 |     vTaskSuspend(serialHandle);
258 |     vTaskSuspend(commandHandle);
259 |     vTaskSuspend(ekgProcessHandle);
260 |     /* Start the scheduler. */
261 |     vTaskStartScheduler();
262 |
263 |     /* Will only get here if there was insufficient memory to create the
idle task.

```

```

264         (Created by vTaskStartScheduler())*/
265     return 0;
266 }
267
268 /*
269  three dummy tasks
270 */
271
272 void measure(void *vParameters)
273 {
274     while(1)
275     {
276         measureTask.runTaskFunction(measureTask.taskDataPtr);
277         vTaskDelay(MINOR_CYCLE * MAJOR_CYCLE);
278     }
279 }
280
281 void compute(void *vParameters)
282 {
283     while(1)
284     {
285         computeTask.runTaskFunction(computeTask.taskDataPtr);
286         vTaskDelay(MINOR_CYCLE);
287     }
288 }
289
290 void keyPad(void *vParameters)
291 {
292     while(1)
293     {
294         keyPadTask.runTaskFunction(keyPadTask.taskDataPtr);
295         vTaskDelay(MINOR_CYCLE);
296     }
297 }
298
299 void display(void *vParameters)
300 {
301     while(1)
302     {
303         displayTask.runTaskFunction(displayTask.taskDataPtr);
304         vTaskDelay(MINOR_CYCLE);
305     }
306 }
307
308 void warning(void *vParameters)
309 {
310     while(1)
311     {
312         warningTask.runTaskFunction(warningTask.taskDataPtr);
313         vTaskDelay(MINOR_CYCLE);
314     }
315 }
316

```

```

317 void serial(void *vParameters)
318 {
319     while(1)
320     {
321         serialTask.runTaskFunction( serialTask.taskDataPtr );
322         vTaskDelay(MINOR_CYCLE);
323     }
324 }
325
326 void status(void *vParameters)
327 {
328     while(1)
329     {
330         statusTask.runTaskFunction( statusTask.taskDataPtr );
331         vTaskDelay(MINOR_CYCLE);
332     }
333 }
334
335 void ekgCapture(void *vParameters)
336 {
337     while(1)
338     {
339         ekgCaptureTask.runTaskFunction( ekgCaptureTask.taskDataPtr );
340         vTaskDelay(MINOR_CYCLE);
341     }
342 }
343
344 void ekgProcess(void *vParameters)
345 {
346     while(1)
347     {
348         ekgProcessTask.runTaskFunction( ekgProcessTask.taskDataPtr );
349         vTaskDelay(MINOR_CYCLE);
350     }
351 }
352
353 void command(void *vParameters)
354 {
355     while(1)
356     {
357         commandTask.runTaskFunction( commandTask.taskDataPtr );
358         vTaskDelay(MINOR_CYCLE);
359     }
360 }
361
362 /*-----*/
363
364 void vApplicationStackOverflowHook( xTaskHandle *pxTask, signed portCHAR *
pcTaskName )
365 {
366     ( void ) pxTask;
367     ( void ) pcTaskName;
368

```



```

369     while( 1 );
370 }
371
372 /*-----*/
373
374 void prvSetupHardware( void )
375 {
376     /*
377      If running on Rev A2 silicon , turn the LDO voltage up to 2.75V.  This
378      is
379      a workaround to allow the PLL to operate reliably.
380     */
381     if( DEVICE_IS_REVA2 )
382     {
383         SysCtlLDOSet( SYSCTL_LDO_2_75V );
384     }
385
386     // Set the clocking to run from the PLL at 50 MHz
387
388     SysCtlClockSet( SYSCTL_SYSDIV_4 | SYSCTL_USE_PLL | SYSCTL_OSC_MAIN |
389 SYSCTL_XTAL_8MHZ );
390
391     /*
392      Enable Port F for Ethernet LEDs
393      LED0      Bit 3      Output
394      LED1      Bit 2      Output
395     */
396
397     SysCtlPeripheralEnable( SYSCTL_PERIPH_GPIOF );
398     GPIODirModeSet( GPIO_PORTF_BASE, ( GPIO_PIN_2 | GPIO_PIN_3 ),
399 GPIO_DIR_MODE_HW );
400     GPIOPadConfigSet( GPIO_PORTF_BASE, ( GPIO_PIN_2 | GPIO_PIN_3 ),
401 GPIO_STRENGTH_2MA, GPIO_PIN_TYPE_STD );
402 }
403
404 /*-----*/
405
406 void vApplicationTickHook( void )
407 {
408     // static xOLEDMMessage xMessage = { "PASS" };
409     static unsigned portLONG ulTicksSinceLastDisplay = 0;
410     portBASE_TYPE xHigherPriorityTaskWoken = pdFALSE;
411
412     /*
413      Called from every tick interrupt.  Have enough ticks passed to make it
414      time to perform our health status check again?
415     */
416
417     ulTicksSinceLastDisplay++;
418     if( ulTicksSinceLastDisplay >= mainCHECK_DELAY )

```

```

418     {
419         ulTicksSinceLastDisplay = 0;
420     }
421 }
422 }

```

## D.2 Global Data

../code/globals.h

```

1  /*
2  * globals.h
3  * Author(s): Jonathan Ellington , Patrick Ma
4  * 1/28/2014
5  *
6  * Defines global data for tasks to access
7  * MUST be initialized before using
8  */
9
10 #include <String.h>
11 #include "inc/hw_types.h"
12 #include "CircularBuffer.h"
13 #include "stdint.h"
14 #include "inc/hw_memmap.h"
15 #include "driverlib/gpio.h"
16 #include "driverlib/sysctl.h"
17 #include <string.h>
18
19 #define TEMP_RAW_INIT 80          // initial 80
20 #define SYS_RAW_INIT 80          // initial 50
21 #define DIA_RAW_INIT 50          // initial 50
22 #define PULSE_RAW_INIT 30        // initial 30
23
24 #define TEMP_CORR_INIT 0.0
25 #define SYS_CORR_INIT 0.0
26 #define DIA_CORR_INIT 0.0
27 #define PULSE_CORR_INIT 0.0
28 #define EKG_FREQ_RLT 0
29
30 #define BATT_INIT 200
31
32 #define NUM_EKG_SAMPLES 256
33 #define SAMPLE_FREQ 8000 // # sample frequency to get a good measure of <
34                          3750 Hz
35 #define COMMAND_LENGTH 10 // length of command string
36 #define RESPONSE_LENGTH 600 // length of response string
37
38 typedef struct global_data {
39     CircularBuffer temperatureRaw;
40     CircularBuffer systolicPressRaw;
41     CircularBuffer diastolicPressRaw;

```

```

42 CircularBuffer pulseRateRaw;
43 int ekgRaw[NUM_EKG_SAMPLES];
44 int ekgTemp[NUM_EKG_SAMPLES];
45
46 CircularBuffer temperatureCorrected;
47 CircularBuffer systolicPressCorrected;
48 CircularBuffer diastolicPressCorrected;
49 CircularBuffer pulseRateCorrected;
50 CircularBuffer ekgFrequencyResult;
51
52 unsigned short batteryState;
53 unsigned short mode;
54 unsigned short measurementSelection;
55 tBoolean measurementComplete;
56 tBoolean ekgCaptureDone;
57 tBoolean ekgProcessDone;
58 tBoolean responseReady;
59 tBoolean alarmAcknowledge;
60 tBoolean select;
61 unsigned short scroll;
62
63 char commandStr[COMMAND_LENGTH];
64 char responseStr[RESPONSE_LENGTH];
65 } GlobalData;
66
67 extern GlobalData global;
68
69 // initializes the global variables for use by system
70 void initializeGlobalData();
71
72 // allows use of pin47 for debug, toggles the pin hi or lo alternatively
73 void debugPin47();

```

../code/globals.c

```

1 /*
2  * globals.c
3  * Author(s): Jonathan Ellington, Patrick Ma
4  * 1/28/2014
5  *
6  * Defines global data for tasks to access
7  */
8 #include "globals.h"
9 #include "utils/ustdlib.h"
10
11 GlobalData global;
12
13 // The arrays to be wrapped in a
14 // circular buffer
15 static int temperatureRawArr[8];
16 static int systolicPressRawArr[8];
17 static int diastolicPressRawArr[8];
18 static int pulseRateRawArr[8];

```

```

19
20 static float temperatureCorrectedArr[8];
21 static float systolicPressCorrectedArr[8];
22 static float diastolicPressCorrectedArr[8];
23 static float pulseRateCorrectedArr[8];
24 static int ekgFrequencyResultArr[16];
25
26 static signed int ekgRaw[NUM_EKG_SAMPLES]; // initialize all the elements
    to 0
27 static signed int ekgTemp[NUM_EKG_SAMPLES];
28
29 void initializeGlobalData() {
30     // Wrap the arrays
31     global.temperatureRaw = cbWrap(temperatureRawArr, sizeof(int), 8);
32     global.systolicPressRaw = cbWrap(systolicPressRawArr, sizeof(int), 8);
33     global.diastolicPressRaw = cbWrap(diastolicPressRawArr, sizeof(int), 8);
34     global.pulseRateRaw = cbWrap(pulseRateRawArr, sizeof(int), 8);
35
36     memset(global.ekgRaw, 0, NUM_EKG_SAMPLES);
37     memset(global.ekgTemp, 0, NUM_EKG_SAMPLES);
38
39     global.temperatureCorrected = cbWrap(temperatureCorrectedArr, sizeof(float
    ), 8);
40     global.systolicPressCorrected = cbWrap(systolicPressCorrectedArr, sizeof(
    float), 8);
41     global.diastolicPressCorrected = cbWrap(diastolicPressCorrectedArr, sizeof
    (float), 8);
42     global.pulseRateCorrected = cbWrap(pulseRateCorrectedArr, sizeof(float),
    8);
43     global.ekgFrequencyResult = cbWrap(ekgFrequencyResultArr, sizeof(int), 16)
    ;
44
45     int tr = TEMP_RAW_INIT;
46     int sr = SYS_RAW_INIT;
47     int dr = DIA_RAW_INIT;
48     int pr = PULSE_RAW_INIT;
49
50     float tc = TEMP_CORR_INIT;
51     float sc = SYS_CORR_INIT;
52     float dc = DIA_CORR_INIT;
53     float pc = PULSE_CORR_INIT;
54     int fr = EKG_FREQ_RLT;
55
56     // Add initial values
57     cbAdd(&(global.temperatureRaw), &tr);
58     cbAdd(&(global.systolicPressRaw), &sr);
59     cbAdd(&(global.diastolicPressRaw), &dr);
60     cbAdd(&(global.pulseRateRaw), &pr);
61
62     cbAdd(&(global.temperatureCorrected), &tc);
63     cbAdd(&(global.systolicPressCorrected), &sc);
64     cbAdd(&(global.diastolicPressCorrected), &dc);
65     cbAdd(&(global.pulseRateCorrected), &pc);

```

```

66  cbAdd(&(global.ekgFrequencyResult), &fr);
67
68  // Set normal variables
69  global.batteryState = 200;
70  global.mode = 1;
71  global.measurementSelection = 0;
72  global.alarmAcknowledge = false;
73  global.select = false;
74  global.scroll = 0;
75
76  memset(&(global.commandStr), NULL, sizeof(char) * COMMAND_LENGTH);
77  memset(&(global.responseStr), NULL, sizeof(char) * RESPONSE_LENGTH);
78
79  global.measurementComplete = false;
80  global.ekgCaptureDone = false;
81  global.ekgProcessDone = false;
82  global.responseReady = false;
83 }
84
85 // debug tool
86 void debugPin47() {
87     SysCtlPeripheralEnable(SYSCTL_PERIPH_GPIOB);           // debug
88     GPIOPinTypeGPIOOutput(GPIO_PORTB_BASE, GPIO_PIN_2);
89     long a = GPIOPinRead(GPIO_PORTB_BASE, GPIO_PIN_2);
90     if (0 == a)
91         GPIOPinWrite(GPIO_PORTB_BASE, GPIO_PIN_2, 0xFF);
92     else
93         GPIOPinWrite(GPIO_PORTB_BASE, GPIO_PIN_2, 0x00);
94 }

```

### D.3 Timebase

../code/timebase.h

```

1  /*
2  * timebase.h
3  * Author(s): Jonathan Ellington
4  * 1/28/2014
5  *
6  * Defines the major and minor cycles the system runs on
7  */
8
9  #include "hw_types.h"
10
11 #define MINOR_CYCLE 250           // minor cycle, in milliseconds
12 #define MAJOR_CYCLE 4             // major cycle, in number of minor cycles
13 #define PULSE_CYCLE 8            // pulse rate measurement time, in number of
14                                 // minor cycles
15 #define IS_PULSE_CYCLE (minor_cycle_ctr % (MAJOR_CYCLE * PULSE_CYCLE) == 0)
16 #define IS_MAJOR_CYCLE (minor_cycle_ctr % MAJOR_CYCLE == 0)
17
18 extern unsigned int minor_cycle_ctr; // counts number of minor cycles

```

```

19
20 // returns whether or not length minor cycles have happened since
21 // start_time
22 tBoolean timeHasPassed(int start_time , int length);
23
24 void initializeTimebase(void);
25 void TimerAIntHandler(void);

```

#### D.4 Modified lmi\_fs.c

../code/lmi\_fs.c

```

1 //
  *****
2 //
3 // lmi_fs.c – File System Processing for enet_io application.
4 //
5 // Copyright (c) 2007–2012 Texas Instruments Incorporated. All rights
  reserved.
6 // Software License Agreement
7 //
8 // Texas Instruments (TI) is supplying this software for use solely and
9 // exclusively on TI’s microcontroller products. The software is owned by
10 // TI and/or its suppliers, and is protected under applicable copyright
11 // laws. You may not combine this software with "viral" open-source
12 // software in order to form a larger program.
13 //
14 // THIS SOFTWARE IS PROVIDED "AS IS" AND WITH ALL FAULTS.
15 // NO WARRANTIES, WHETHER EXPRESS, IMPLIED OR STATUTORY, INCLUDING, BUT
16 // NOT LIMITED TO, IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR
17 // A PARTICULAR PURPOSE APPLY TO THIS SOFTWARE. TI SHALL NOT, UNDER ANY
18 // CIRCUMSTANCES, BE LIABLE FOR SPECIAL, INCIDENTAL, OR CONSEQUENTIAL
19 // DAMAGES, FOR ANY REASON WHATSOEVER.
20 //
21 // This is part of revision 9107 of the EK-LM3S8962 Firmware Package.
22 //
23 //
  *****
24
25 #include <string.h>
26 #include "inc/hw_types.h"
27 #include "utils/lwiplib.h"
28 #include "utils/ustdlib.h"
29 #include "httpserver_raw/fs.h"
30 #include "httpserver_raw/fsdata.h"
31
32 #include "FreeRTOS.h"
33 #include "task.h"
34 #include "globals.h"
35

```

```

36 extern xTaskHandle commandHandle;
37 // #include "io.h"
38
39 //
    *****
40 //
41 // Include the file system data for this application. This file is
    generated
42 // by the makefsfile utility , using the following command (all on one line):
43 //
44 //     makefsfile -i fs -o io_fsdata.h -r -h
45 //
46 // If any changes are made to the static content of the web pages served by
    the
47 // application , this script must be used to regenerate io_fsdata.h in order
48 // for those changes to be picked up by the web server.
49 //
50 //
    *****
51 #include "io_fsdata.h"
52
53 //
    *****
54 //
55 // Global Settings for demo page content.
56 //
57 //
    *****
58 static char g_cSampleTextBuffer[10] = {0};
59
60 //
    *****
61 //
62 // Open a file and return a handle to the file , if found. Otherwise ,
63 // return NULL. This function also looks for special filenames used to
64 // provide specific status information or to control various subsystems.
65 // These filenames are used by the JavaScript on the "IO Control Demo 1"
66 // example web page.
67 //
68 //
    *****
69 struct fs_file *
70 fs_open(char *name)
71 {
72     char *data;
73     int i;
74     const struct fsdata_file *ptTree;

```

```

75 struct fs_file *ptFile = NULL;
76
77 //
78 // Allocate memory for the file system structure.
79 //
80 ptFile = mem_malloc(sizeof(struct fs_file));
81 if(NULL == ptFile)
82 {
83     return(NULL);
84 }
85
86 //
87 // Process command request
88 //
89 if (strcmp(name, "/cgi-bin/send_command/value=", 28) == 0)
90 {
91     // Get Command String
92     data = name;
93     data += 28;
94     i = 0;
95     do
96     {
97         switch(data[i])
98         {
99             case '+':
100                 global.commandStr[i] = ' ';
101                 break;
102             default:
103                 global.commandStr[i] = data[i];
104                 break;
105         }
106         if(global.commandStr[i] == 0)
107         {
108             break;
109         }
110         i++;
111     }while(i < sizeof(global.commandStr));
112
113     // Setup the file structure to return whatever.
114     ptFile->data = NULL;
115     ptFile->len = 0;
116     ptFile->index = 0;
117     ptFile->pextension = NULL;
118
119     vTaskResume(commandHandle);
120     return(ptFile);
121 }
122
123 //
124 // Process command response
125 //
126 if (strcmp(name, "/cgi-bin/receive_command", 24) == 0)
127 {

```



```

128 if (global.responseReady) {
129
130     char *buf = (char *) mem_malloc(sizeof(global.responseStr));
131     memcpy(buf, global.responseStr, sizeof(global.responseStr));
132
133     // Setup the file structure to return whatever.
134     ptFile->data = buf;
135     ptFile->len = strlen(buf);
136     ptFile->index = 0;
137     ptFile->pextension = NULL;
138
139     //
140     // Return the file system pointer.
141     //
142     return(ptFile);
143 }
144
145 ptFile->data = NULL;
146 ptFile->len = 0;
147 ptFile->index = 0;
148 ptFile->pextension = NULL;
149
150 return ptFile;
151 }
152
153 //
154 // If I can't find it there, look in the rest of the main file system
155 //
156 else
157 {
158     //
159     // Initialize the file system tree pointer to the root of the linked
160     // list.
161     //
162     ptTree = FS_ROOT;
163
164     //
165     // Begin processing the linked list, looking for the requested file
166     // name.
167     //
168     while(NULL != ptTree)
169     {
170         //
171         // Compare the requested file "name" to the file name in the
172         // current node.
173         //
174         if(strncmp(name, (char *)ptTree->name, ptTree->len) == 0)
175         {
176             //
177             // Fill in the data pointer and length values from the
178             // linked list node.
179             //
180             ptFile->data = (char *)ptTree->data;

```

```

179         ptFile->len = ptTree->len;
180
181         //
182         // For now, we setup the read index to the end of the file ,
183         // indicating that all data has been read.
184         //
185         ptFile->index = ptTree->len;
186
187         //
188         // We are not using any file system extensions in this
189         // application , so set the pointer to NULL.
190         //
191         ptFile->pextension = NULL;
192
193         //
194         // Exit the loop and return the file system pointer.
195         //
196         break;
197     }
198
199     //
200     // If we get here, we did not find the file at this node of the
linked
201     // list.  Get the next element in the list.
202     //
203     ptTree = ptTree->next;
204 }
205 }
206
207 //
208 // If we didn't find the file , ptTee will be NULL.  Make sure we
209 // return a NULL pointer if this happens.
210 //
211 if(NULL == ptTree)
212 {
213     mem_free(ptFile);
214     ptFile = NULL;
215 }
216
217 //
218 // Return the file system pointer.
219 //
220 return(ptFile);
221 }
222
223 //
*****
224 //
225 // Close an opened file designated by the handle.
226 //
227 //
*****

```

```

228 void
229 fs_close(struct fs_file *file)
230 {
231     //
232     // Free the main file system object.
233     //
234     mem_free(file);
235 }
236
237 //
238 // *****
239 // Read the next chunk of data from the file. Return the count of data
240 // that was read. Return 0 if no data is currently available. Return
241 // a -1 if at the end of file.
242 //
243 // *****
244
245 int
246 fs_read(struct fs_file *file, char *buffer, int count)
247 {
248     int iAvailable;
249
250     //
251     // Check to see if a command (pextension = 1).
252     //
253     if(file->pextension == (void *)1)
254     {
255         //
256         // Nothing to do for this file type.
257         //
258         file->pextension = NULL;
259         return(-1);
260     }
261
262     //
263     // Check to see if more data is available.
264     //
265     if(file->len == file->index)
266     {
267         //
268         // There is no remaining data. Return a -1 for EOF indication.
269         //
270         return(-1);
271     }
272
273     //
274     // Determine how much data we can copy. The minimum of the 'count'
275     // parameter or the available data in the file system buffer.
276     //

```

```

276     iAvailable = file->len - file->index;
277     if(iAvailable > count)
278     {
279         iAvailable = count;
280     }
281
282     //
283     // Copy the data.
284     //
285     memcpy(buffer, file->data + iAvailable, iAvailable);
286     file->index += iAvailable;
287
288     //
289     // Return the count of data that we copied.
290     //
291     return(iAvailable);
292 }

```

## D.5 Circular Buffer

../code/CircularBuffer.h

```

1  /* CircularBuffer.h
2   * Jonathan Ellington
3   * 2/7/14
4   */
5
6  /* A circular buffer implementation.
7   * Meant to work without any dynamically allocated memory by wrap()ing
8   * user defined arrays.
9   *
10  * NOTE: This implementation keeps NO type information. This means the
11  *       USER IS RESPONSIBLE FOR KEEPING TRACK OF THE TYPE OF ARRAY THAT
12  *       HAS BEEN WRAPPED!
13  *
14  * NOTE: The implementation here should also be hidden, but this is
15  *       troublesome without dynamic memory. Do not rely on any of these
16  *       elements!
17  */
18
19 #ifndef _CIRCULAR_BUFFER_H
20 #define _CIRCULAR_BUFFER_H
21
22 typedef struct _circBuf {
23     void *array;
24     int sizeElm;
25     int nElm;
26     int currElm;
27 } CircularBuffer;
28
29 /* Wrap an array in a circular buffer.
30  *

```

```

31  * @param arr      the array to wrap
32  * @param sizeElem the size of each element, in bytes
33  * @param nElem    the number of elements in the array
34  *
35  * This function expects the array is freshly created. It will overwrite
36  * array contents on adds.
37  */
38  CircularBuffer cbWrap(void *arr, int sizeElem, int nElem);
39
40  /* Returns a pointer to the current element in the circular buffer
41   * Be sure not to clobber this value! */
42  void *cbGet(CircularBuffer *cb);
43
44  /* Returns a pointer to the wrapped array, with the oldest element
45   * at the end and the newest element at the beginning */
46  void *cbGetArray(CircularBuffer *cb);
47
48  /* Adds elem to cb */
49  void cbAdd(CircularBuffer *cb, void *elem);
50
51  #endif // _CIRCULAR_BUFFER_H

```

../code/CircularBuffer.c

```

1  /* CircularBuffer.c
2   * Jonathan Ellington
3   * 2/7/14
4   */
5
6  #include "CircularBuffer.h"
7  #include <stdio.h>
8
9  CircularBuffer cbWrap(void *arr, int se, int ne) {
10     CircularBuffer cb;
11     cb.array = arr;
12     cb.sizeElm = se;
13     cb.nElm = ne;
14     cb.currElm = 0;
15
16     return cb;
17 }
18
19 /* Returns the current element in the circular buffer */
20 void *cbGet(CircularBuffer *cb) {
21     int sizeElm = cb->sizeElm;
22     int index = cb->currElm;
23
24     unsigned char *bytePtr = (unsigned char *)cb->array;
25     bytePtr += sizeElm * index;
26
27     return (void *) bytePtr;
28 }
29

```

```

30 /* Returns a pointer to the wrapped array, with the oldest element
31 * at the end and the newest element at the beginning */
32 void *cbGetArray(CircularBuffer *cb) {
33     return NULL;
34 }
35
36 /* Adds elm to cb */
37 void cbAdd(CircularBuffer *cb, void *elm) {
38     cb->currElm = (cb->currElm + 1) % cb->nElm;
39
40     int sizeElm = cb->sizeElm;
41     int index = cb->currElm;
42
43     // copy sizeElm bytes from elm into the right spot
44     unsigned char *elmPtr = (unsigned char *) elm;
45     unsigned char *arrElmPtr = (unsigned char *) cb->array;
46     arrElmPtr += sizeElm * index;
47
48     for (int i = 0; i < sizeElm; i++) {
49         arrElmPtr[i] = elmPtr[i];
50     }
51 }

```

## D.6 Warm up File (For Interrupt Initialization)

../code/startup\_ewarm.c

```

1 //
2 //
3 // startup_ewarm.c – Boot code for Stellaris.
4 //
5 // Copyright (c) 2006–2007 Luminary Micro, Inc. All rights reserved.
6 //
7 // Software License Agreement
8 //
9 // Luminary Micro, Inc. (LMI) is supplying this software for use solely and
10 // exclusively on LMI's microcontroller products.
11 //
12 // The software is owned by LMI and/or its suppliers, and is protected under
13 // applicable copyright laws. All rights are reserved. Any use in
14 // violation
15 // of the foregoing restrictions may subject the user to criminal sanctions
16 // under applicable laws, as well as to civil liability for the breach of
17 // the
18 // terms and conditions of this license.
19 //
20 // THIS SOFTWARE IS PROVIDED "AS IS". NO WARRANTIES, WHETHER EXPRESS,
21 // IMPLIED
22 // OR STATUTORY, INCLUDING, BUT NOT LIMITED TO, IMPLIED WARRANTIES OF

```

```

20 // MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE APPLY TO THIS
    // SOFTWARE.
21 // LMI SHALL NOT, IN ANY CIRCUMSTANCES, BE LIABLE FOR SPECIAL, INCIDENTAL,
    // OR
22 // CONSEQUENTIAL DAMAGES, FOR ANY REASON WHATSOEVER.
23 //
24 // This is part of revision 100 of the Stellaris Ethernet
25 // Applications Library.
26 //
27 //
    *****
28
29 //
    *****
30 //
31 // Enable the IAR extensions for this source file.
32 //
33 //
    *****
34 #pragma language=extended
35
36 //
    *****
37 //
38 // Forward declaration of the default fault handlers.
39 //
40 //
    *****
41 static void NmiSR(void);
42 static void FaultISR(void);
43 static void IntDefaultHandler(void);
44 extern void lwIPEthernetIntHandler(void);
45 extern void timer1IntHandler(void);
46
47 //
    *****
48 //
49 // External declaration for the interrupt handler used by the application.
50 //
51 //
    *****
52
53
54 //
    *****

```

```

55 //
56 // The entry point for the application.
57 //
58 //
    *****

59 extern void __iar_program_start(void);
60 extern void xPortPendSVHandler(void);
61 extern void xPortSysTickHandler(void);
62 extern void vPortSVCHandler(void);
63 extern void vT2InterruptHandler( void );
64 extern void vT3InterruptHandler( void );
65 extern void vEMAC_ISR( void );
66 //extern void Timer0IntHandler( void ); not sure where this is from
67 extern void TimerAIntHandler( void ); // timebase.h
68 extern void ADC0IntHandler(void); // ADC0 — for EKG sensor
69
70 //
    *****

71 //
72 // Reserve space for the system stack.
73 //
74 //
    *****

75 #ifndef STACK_SIZE
76 #define STACK_SIZE 1024
77 #endif
78 static unsigned long pulStack[STACK_SIZE] @ ".noinit";
79
80 //
    *****

81 //
82 // A union that describes the entries of the vector table. The union is
    needed
83 // since the first entry is the stack pointer and the remainder are function
84 // pointers.
85 //
86 //
    *****

87 typedef union
88 {
89     void (*pfnHandler)(void);
90     unsigned long ulPtr;
91 }
92 uVectorEntry;
93
94 //
    *****

```



```

95 //
96 // The minimal vector table for a Cortex-M3. Note that the proper
   constructs
97 // must be placed on this to ensure that it ends up at physical address
98 // 0x0000.0000.
99 //
100 //
   *****

101 __root const uVectorEntry __vector_table[] @ ".intvec" =
102 {
103     { .ulPtr = (unsigned long)pulStack + sizeof(pulStack) },
104                                     // The initial stack pointer
105     __iar_program_start,           // The reset handler
106     NmiSR,                         // The NMI handler
107     FaultISR,                     // The hard fault handler
108     IntDefaultHandler,            // The MPU fault handler
109     IntDefaultHandler,            // The bus fault handler
110     IntDefaultHandler,            // The usage fault handler
111     0,                            // Reserved
112     0,                            // Reserved
113     0,                            // Reserved
114     0,                            // Reserved
115     vPortSVCHandler,              // SVC call handler
116     IntDefaultHandler,            // Debug monitor handler
117     0,                            // Reserved
118     xPortPendSVHandler,           // The PendSV handler
119     xPortSysTickHandler,          // The SysTick handler
120     IntDefaultHandler,            // GPIO Port A
121     IntDefaultHandler,            // GPIO Port B
122     IntDefaultHandler,            // GPIO Port C
123     IntDefaultHandler,            // GPIO Port D
124     IntDefaultHandler,            // GPIO Port E
125     IntDefaultHandler,            // UART0 Rx and Tx
126     IntDefaultHandler,            // UART1 Rx and Tx
127     IntDefaultHandler,            // SSI Rx and Tx
128     IntDefaultHandler,            // I2C Master and Slave
129     IntDefaultHandler,            // PWM Fault
130     IntDefaultHandler,            // PWM Generator 0
131     IntDefaultHandler,            // PWM Generator 1
132     IntDefaultHandler,            // PWM Generator 2
133     IntDefaultHandler,            // Quadrature Encoder
134     IntDefaultHandler,            // ADC Sequence 0
135     ADC0IntHandler,               // ADC Sequence 1
136     IntDefaultHandler,            // ADC Sequence 2
137     IntDefaultHandler,            // ADC Sequence 3
138     IntDefaultHandler,            // Watchdog timer
139     TimerAIntHandler,             // Timer 0 subtimer A
140     IntDefaultHandler,            // Timer 0 subtimer B
141     IntDefaultHandler,            // Timer 1 subtimer A
142     IntDefaultHandler,            // Timer 1 subtimer B
143     IntDefaultHandler,            // Timer 2 subtimer A
144     IntDefaultHandler,            // Timer 2 subtimer B

```

```

145     IntDefaultHandler ,           // Analog Comparator 0
146     IntDefaultHandler ,           // Analog Comparator 1
147     IntDefaultHandler ,           // Analog Comparator 2
148     IntDefaultHandler ,           // System Control (PLL, OSC, BO)
149     IntDefaultHandler ,           // FLASH Control
150     IntDefaultHandler ,           // GPIO Port F
151     IntDefaultHandler ,           // GPIO Port G
152     IntDefaultHandler ,           // GPIO Port H
153     IntDefaultHandler ,           // UART2 Rx and Tx
154     IntDefaultHandler ,           // SSI1 Rx and Tx
155     IntDefaultHandler ,           // Timer 3 subtimer A
156     IntDefaultHandler ,           // Timer 3 subtimer B
157     IntDefaultHandler ,           // I2C1 Master and Slave
158     IntDefaultHandler ,           // Quadrature Encoder 1
159     IntDefaultHandler ,           // CAN0
160     IntDefaultHandler ,           // CAN1
161     IntDefaultHandler ,           // CAN2
162     lwIPEthernetIntHandler ,       // Ethernet
163     IntDefaultHandler ,           // Hibernate
164     IntDefaultHandler ,           // USB0
165     IntDefaultHandler ,           // PWM Generator 3
166     IntDefaultHandler ,           // uDMA Software Transfer
167     IntDefaultHandler ,           // uDMA Error
168 };
169
170
171 //
172 // *****
173 // This is the code that gets called when the processor receives a NMI.
174 // This
175 // simply enters an infinite loop, preserving the system state for
176 // examination
177 // by a debugger.
178 //
179 // *****
180
181 static void
182 NmiSR( void )
183 {
184     //
185     // Enter an infinite loop.
186     //
187     while(1)
188     {
189     }
190 }
191 //
192 // *****

```

```

190 //
191 // This is the code that gets called when the processor receives a fault
192 // interrupt. This simply enters an infinite loop, preserving the system
    state
193 // for examination by a debugger.
194 //
195 //
    *****

196 static void
197 FaultISR(void)
198 {
199     //
200     // Enter an infinite loop.
201     //
202     while(1)
203     {
204     }
205 }
206 //
207 //
    *****

208 //
209 // This is the code that gets called when the processor receives an
    unexpected
210 // interrupt. This simply enters an infinite loop, preserving the system
    state
211 // for examination by a debugger.
212 //
213 //
    *****

214 static void
215 IntDefaultHandler(void)
216 {
217     //
218     // Go into an infinite loop.
219     //
220     while(1)
221     {
222     }
223 }

```

## D.7 Tasks

### D.7.1 Task Control Blocks

../code/tcb.h

```

1 /*
2  * task.h

```

```

3  * Author(s): Jonathan Ellington
4  * 1/28/2014
5  *
6  * Defines the interface for a task
7  */
8
9  #ifndef _TASK_H
10 #define _TASK_H
11
12 typedef struct tcb_struct {
13     void (*runTaskFunction) (void*);
14     void *taskDataPtr;
15 } TCB;
16
17 #endif // _TASK_H

```

### D.7.2 Startup Task

../code/startup.h

```

1  /*
2  * startup.h
3  * author: Patrick ma
4  * Date: 2/13/2014
5  *
6  * initializes the system and system variables
7  */
8
9  extern tBoolean runSchedule;
10
11 // Interrupt handler for hw timer
12 void SysTickIntHandler();
13
14 // configure the hardware timer
15 void initializeHWCounter();
16
17 // initialize system variables and hardware timer
18 void startup();
19
20 void DisplayIPAddress();
21
22 void lwIPHostTimerHandler();

```

../code/startup.c

```

1  /*
2  * startup.c
3  * author: Patrick ma
4  * Date: 2/13/2014
5  *
6  * initializes the system and system variables
7  */

```

```

8
9 #include <String.h>
10 #include "timebase.h"
11 #include "globals.h"
12 #include "driverlib/systick.h" // for systick interrupt (minor cycle)
13 #include "driverlib/sysctl.h"
14 #include "driverlib/interrupt.h"
15 #include "startup.h"
16 #include "schedule.h"
17 #include "inc/hw_ints.h"
18 #include "inc/hw_memmap.h"
19 #include "inc/hw_nvic.h"
20 #include "inc/hw_types.h"
21 #include "driverlib/ethernet.h"
22 #include "driverlib/flash.h"
23 #include "driverlib/gpio.h"
24 #include "utils/locator.h"
25 #include "utils/lwiplib.h"
26 #include "utils/uartstdio.h"
27 #include "utils/ustdlib.h"
28 #include "httpserver_raw/httpd.h"
29 #include "drivers/rit128x96x4.h"
30 #include "driverlib/timer.h"
31
32
33
34 //
35 // *****
36 // Defines for setting up the system clock.
37 //
38 //
39 // *****
40
41 #define SYSTICKHZ          100
42 #define SYSTICKMS          (1000 / SYSTICKHZ)
43 #define SYSTICKUS          (1000000 / SYSTICKHZ)
44 #define SYSTICKNS          (1000000000 / SYSTICKHZ)
45
46 // *****
47
48 //
49 // A set of flags. The flag bits are defined as follows:
50 //
51 // 0 -> An indicator that a SysTick interrupt has occurred.
52 //
53 // *****
54
55 #define FLAG_SYSTICK          0
56 static volatile unsigned long g_ulFlags;

```

```

53 |
54 | //
55 | // *****
56 | // External Application references.
57 | //
58 | // *****
59 | extern void httpd_init(void);
60 |
61 |
62 | //
63 | // *****
64 | // Timeout for DHCP address request (in seconds).
65 | //
66 | // *****
67 | #ifndef DHCP_EXPIRE_TIMER_SECS
68 | #define DHCP_EXPIRE_TIMER_SECS 45
69 | #endif
70 |
71 | void timer1IntHandler (void) {
72 |     TimerIntClear(TIMER1_BASE, TIMER_TIMA_TIMEOUT);
73 |     //
74 |     // Indicate that a SysTick interrupt has occurred.
75 |     //
76 |     HWREGBITW(&g_uIFlags, FLAG_SYSTICK) = 1;
77 |     //
78 |     // Call the lwIP timer handler.
79 |     //
80 |     lwIPTimer(SYSTICKMS);
81 | }
82 |
83 |
84 |
85 |
86 | /*
87 | * Initializes hw counter and system state variables
88 | */
89 | void startup() {
90 |     SysCtlClockSet(SYSCTL_SYSDIV_4 | SYSCTL_USE_PLL | SYSCTL_OSC_MAIN |
91 |     SYSCTL_XTAL_8MHZ);
92 |
93 |
94 |
95 |     RIT128x96x4Init(1000000);
96 |
97 |     // Initialize global data

```

```

98  initializeGlobalData();    // from globals.h
99  initializeTimebase();
100
101
102
103  //////////////////////////////////////////////////
104  // Web server initialization stuff //
105  //////////////////////////////////////////////////
106  unsigned long ulUser0, ulUser1;
107  unsigned char pucMACArray[8];
108
109  // Enable and Reset the Ethernet Controller.
110  SysCtlPeripheralEnable(SYSCTL_PERIPH_ETH);
111  SysCtlPeripheralReset(SYSCTL_PERIPH_ETH);
112
113  // Enable Port F for Ethernet LEDs.
114  // LED0      Bit 3   Output
115  // LED1      Bit 2   Output
116  SysCtlPeripheralEnable(SYSCTL_PERIPH_GPIOF);
117  GPIOPinTypeEthernetLED(GPIO_PORTF_BASE, GPIO_PIN_2 | GPIO_PIN_3);
118
119  // Configure timer for a periodic interrupt.
120  SysCtlPeripheralEnable(SYSCTL_PERIPH_TIMER1);
121  TimerDisable(TIMER1_BASE, TIMER_BOTH);
122  TimerConfigure(TIMER1_BASE, TIMER_CFG_32_BIT_PER );
123  TimerLoadSet(TIMER1_BASE, TIMER_A, (SysCtlClockGet() / SYSTICKHZ));
124  TimerIntRegister(TIMER1_BASE, TIMER_A, timer1IntHandler);
125  TimerIntEnable(TIMER1_BASE, TIMER_TIMA_TIMEOUT);
126  TimerEnable(TIMER1_BASE, TIMER_A);
127
128  IntMasterEnable();
129
130  FlashUserGet(&ulUser0, &ulUser1);
131  if((ulUser0 == 0xffffffff) || (ulUser1 == 0xffffffff))
132  {
133      // We should never get here. This is an error if the MAC address
134      // has not been programmed into the device. Exit the program.
135      RIT128x96x4StringDraw("MAC Address", 0, 76, 15);
136      RIT128x96x4StringDraw("Not Programmed!", 0, 86, 15);
137      while(1);
138  }
139
140  pucMACArray[0] = ((ulUser0 >> 0) & 0xff);
141  pucMACArray[1] = ((ulUser0 >> 8) & 0xff);
142  pucMACArray[2] = ((ulUser0 >> 16) & 0xff);
143  pucMACArray[3] = ((ulUser1 >> 0) & 0xff);
144  pucMACArray[4] = ((ulUser1 >> 8) & 0xff);
145  pucMACArray[5] = ((ulUser1 >> 16) & 0xff);
146
147  // Initialize the lwIP library, using DHCP.
148  lwIPInit(pucMACArray, 0, 0, 0, IPADDR_USE_DHCP);
149
150  // Setup the device locator service.

```

```

151 LocatorInit();
152 LocatorMACAddrSet(pucMACArray);
153 LocatorAppTitleSet("EK-LM3S8962 enet_io");
154
155 // Initialize a sample httpd server.
156 httpd_init();
157 }
158
159 //
160 // *****
161 // Display an lwIP type IP Address.
162 //
163 //
164 // *****
165 void DisplayIPAddress(unsigned long ipaddr, unsigned long ulCol,
166                      unsigned long ulRow)
167 {
168     char pucBuf[16];
169     unsigned char *pucTemp = (unsigned char *)&ipaddr;
170
171     //
172     // Convert the IP Address into a string.
173     //
174     usprintf(pucBuf, "%d.%d.%d.%d", pucTemp[0], pucTemp[1], pucTemp[2],
175             pucTemp[3]);
176
177     //
178     // Display the string.
179     //
180     RIT128x96x4StringDraw(pucBuf, ulCol, ulRow, 15);
181 }
182 //
183 // *****
184 // Required by lwIP library to support any host-related timer functions.
185 //
186 //
187 // *****
188 void lwIPHostTimerHandler(void)
189 {
190     static unsigned long ulLastIPAddress = 0;
191     unsigned long ulIPAddress;
192
193     ulIPAddress = lwIPLocalIPAddrGet();
194
195     //

```



```

195 // If IP Address has not yet been assigned, update the display
196 // accordingly
197 //
198 if (ulIPAddress == 0)
199 {
200     static int iColumn = 6;
201     //
202     // Update status bar on the display.
203     //
204 //     RIT128x96x4Enable(1000000);
205     if (iColumn < 12)
206     {
207         RIT128x96x4StringDraw(">", 114, 86, 15);
208         RIT128x96x4StringDraw("<", 0, 86, 15);
209         RIT128x96x4StringDraw("*", iColumn, 86, 7);
210     }
211     else
212     {
213         RIT128x96x4StringDraw("*", iColumn - 6, 86, 7);
214     }
215     iColumn += 4;
216     if (iColumn > 114)
217     {
218         iColumn = 6;
219         RIT128x96x4StringDraw(">", 114, 86, 15);
220     }
221     // RIT128x96x4Disable();
222 }
223 //
224 // Check if IP address has changed, and display if it has.
225 //
226 //
227 //
228 else if (ulLastIPAddress != ulIPAddress)
229 {
230     ulLastIPAddress = ulIPAddress;
231 //     RIT128x96x4Enable(1000000);
232 //     RIT128x96x4StringDraw(" ", 0, 16, 15);
233 RIT128x96x4StringDraw(" ", 0, 86, 15);
234 RIT128x96x4StringDraw("IP: ", 0, 86, 15);
235 //     RIT128x96x4StringDraw("MASK: ", 0, 24, 15);
236 //     RIT128x96x4StringDraw("GW: ", 0, 32, 15);
237 DisplayIPAddress(ulIPAddress, 36, 86);
238 //     ulIPAddress = lwIPLocalNetMaskGet();
239 //     DisplayIPAddress(ulIPAddress, 36, 24);
240 //     ulIPAddress = lwIPLocalGWAddrGet();
241 //     DisplayIPAddress(ulIPAddress, 36, 32);
242 //     RIT128x96x4Disable();
243 }
244 }

```

### D.7.3 Measure Task

../code/measure.h

```
1  /*
2  * measure.h
3  * Author(s): Jonathan Ellington
4  * 1/28/2014
5  *
6  * Defines the interface for the measureTask.
7  * initializeMeasureData() should be called before running measureTask()
8  */
9
10 #include "tcb.h"
11
12 /* Points to the TCB for measure */
13 extern TCB measureTask;
```

../code/measure.c

```
1  /*
2  * measure.c
3  * Author(s): Jonathan Ellington
4  * 1/28/2014
5  *
6  * Implements measure.c
7  *
8  * Uses port PA4 for interrupt input for pulse rate
9  */
10
11 #define DEBUG_MEASURE 0
12
13 #include "FreeRTOS.h"
14 #include "task.h"
15
16 #include "CircularBuffer.h"
17 #include "globals.h"
18 #include "timebase.h"
19 #include "measure.h"
20 #include "schedule.h"
21
22 #include "inc/hw_memmap.h"
23 #include "inc/hw_types.h"
24 #include "driverlib/gpio.h"
25 #include "driverlib/sysctl.h"
26 #include "driverlib/interrupt.h"
27 #include "inc/hw_ints.h"
28 #include "driverlib/debug.h"
29 #include "driverlib/adc.h"
30
31 // Used for debug display
32 #if DEBUG_MEASURE
33 #include "drivers/rit128x96x4.h"
34 #include "utils/ustdlib.h"
```

```

35 #include "compute.h"
36 #include "ekgCapture.h"
37 #endif
38
39 // prototype for compiler
40 void measureRunFunction(void *dataptr);
41 void PulseRateISR(void);
42
43 // Internal data structure
44 typedef struct measureData {
45     CircularBuffer *temperatureRaw;
46     CircularBuffer *systolicPressRaw;
47     CircularBuffer *diastolicPressRaw;
48     CircularBuffer *pulseRateRaw;
49     unsigned short *measureSelect;
50     tBoolean *ekgCaptureDone;
51 } MeasureData;
52
53 static int pulseRate = 0;
54 static MeasureData data; // internal data
55 TCB measureTask = {&measureRunFunction, &data}; // task interface
56 extern xTaskHandle ekgCaptureHandle;
57
58 void initializeMeasureTask() {
59     // Load data memory
60     data.temperatureRaw = &(global.temperatureRaw);
61     data.systolicPressRaw = &(global.systolicPressRaw);
62     data.diastolicPressRaw = &(global.diastolicPressRaw);
63     data.pulseRateRaw = &(global.pulseRateRaw);
64     data.measureSelect = &(global.measurementSelection);
65     data.ekgCaptureDone = &(global.ekgCaptureDone);
66
67     //setup for temperature sensor
68     ADCSequenceDisable(ADC0_BASE, 1);
69     ADCSequenceConfigure(ADC0_BASE, 1, ADC_TRIGGER_PROCESSOR, 1);
70     ADCSequenceStepConfigure(ADC0_BASE, 1, 0, ADC_CTL_END | ADC_CTL_TS);
71
72     /* Interrupt setup
73     * Note: using statically registered interrupts, because they're faster
74     *       this means we aren't using the dynamic GPIOPortIntRegister()
75     *       function,
76     *       instead, an entry in the interrupt table is populated with the
77     *       address
78     *       of the interrupt handler (under GPIO Port A) and this is enabled
79     *       with
80     *       IntEnable(INT_GPIOA) */
81
82     SysCtlPeripheralEnable(SYSCTL_PERIPH_GPIOA);
83
84     // Set PA4 as input
85     GPIOPinTypeGPIOInput(GPIO_PORTA_BASE, GPIO_PIN_4);
86
87     // Enable interrupts on PA4

```

```

85  GPIOPinIntEnable(GPIO_PORTA_BASE, GPIO_PIN_4);
86
87  // Set interrupt type
88  GPIOIntTypeSet(GPIO_PORTA_BASE, GPIO_PIN_4, GPIO_RISING_EDGE);
89
90  // Enable interrupts to the processor.
91  IntMasterEnable();
92
93  // Enable the interrupts.
94  IntEnable(INT_GPIOA);
95 }
96
97 void setTemp(CircularBuffer *tbuf) {
98     int temp = 0;
99     ADCSequenceEnable(ADC0_BASE, 1);
100
101     // Trigger the sample sequence.
102     ADCProcessorTrigger(ADC0_BASE, 1);
103
104     // Wait until the sample sequence has completed.
105     while( 1 != ADCSequenceDataGet(ADC0_BASE, 1, &temp));
106
107     ADCSequenceDisable(ADC0_BASE, 1);
108     cbAdd(tbuf, &temp);
109 }
110
111 void setBloodPress(CircularBuffer *spbuf, CircularBuffer *dpbuf) {
112     // This is written to lab spec, with a flag to indicate "complete".
113     // Right now, it does nothing, but I imagine it should probably be a
114     // global
115     // variable to indicate to the compute task that the pressure measurement
116     // is ready, since this measurement takes a nontrivial amount of time
117
118     static unsigned int i = 0;
119
120     static tBoolean sysComplete = false;
121     static tBoolean diaComplete = false;
122
123     int syspress = *(int *)cbGet(spbuf);
124
125     // Restart systolic measurement if diastolic is complete
126
127     if (syspress > 100)
128         sysComplete = true;
129
130     if (! sysComplete) {
131         if (i%2==0) syspress+=3;
132         else syspress--;
133     }
134
135     int diapress = *(int *)cbGet(dpbuf);
136
137     // Restart diastolic measurement if systolic is complete

```

```

137 | if (diapress < 40)
138 |     diaComplete = true;
139 |
140 | if (!diaComplete) {
141 |     if (i%2==0) diapress -=2;
142 |     else diapress++;
143 | }
144 |
145 | if (diaComplete && sysComplete) {
146 |     sysComplete = false;
147 |     diaComplete = false;
148 |     syspress = SYS_RAW_INIT;
149 |     diapress = DIA_RAW_INIT;
150 | }
151 |
152 | cbAdd(spbuf, &syspress);
153 | cbAdd(dpbuf, &diapress);
154 |
155 | i++;
156 | }
157 |
158 | // pulse rate interrupt handler
159 | void PulseRateISR(void) {
160 |     pulseRate++;
161 |
162 |     GPIOPinIntClear(GPIO_PORTA_BASE, GPIO_PIN_4);
163 | }
164 |
165 | void measureRunFunction(void *dataptr) {
166 |     static tBoolean onFirstRun = true;
167 |     static int rate;
168 |     MeasureData *mData = (MeasureData *) dataptr;
169 |
170 |     if (onFirstRun) {
171 |         initializeMeasureTask();
172 |         rate = *(int*) cbGet(mData->pulseRateRaw);
173 |         onFirstRun = false;
174 |     }
175 |
176 |     // capture pulse rate
177 |     if (IS_PULSE_CYCLE) {
178 |         // Divide by two so raw pulse rate matches frequency
179 |         rate = pulseRate/2;
180 |         pulseRate = 0;
181 |     }
182 |
183 |     // only run on major cycle
184 |     short measureSelect = *(mData->measureSelect);
185 |     // if (measureSelect == 0 || measureSelect == 1)
186 |     // {
187 |     setTemp(mData->temperatureRaw);
188 |     // }
189 |     // if (measureSelect == 0 || measureSelect == 2)

```

```

190 // {
191     setBloodPress(mData->systolicPressRaw, mData->diastolicPressRaw);
192 // }
193 // if (measureSelect == 0 || measureSelect == 3)
194 // {
195     int prev = *(int*) cbGet(mData->pulseRateRaw);
196
197     // Only save if +- 15%
198     if (rate < prev*0.85 || rate > prev*1.15) {
199         cbAdd(mData->pulseRateRaw, (void *)&rate);
200     }
201 // }
202 // if (measureSelect == 0 || measureSelect == 4)
203 // {
204     *(mData->ekgCaptureDone) = false;
205     vTaskResume(ekgCaptureHandle);
206 #if DEBUG_MEASURE
207     RIT128x96x4StringDraw("ekgCapture go!", 0, 50, 15);
208 #endif
209 // }
210 // else
211 // {
212     vTaskSuspend(ekgCaptureHandle);
213 #if DEBUG_MEASURE
214     RIT128x96x4StringDraw("no ekg!", 0, 50, 15);
215 #endif
216 // }
217
218 #if DEBUG_MEASURE
219     char num[30];
220     int temp = *(int *)cbGet(mData->temperatureRaw);
221     int sys = *(int *)cbGet(mData->systolicPressRaw);
222     int dia = *(int *)cbGet(mData->diastolicPressRaw);
223     int pulse = *(int *)cbGet(mData->pulseRateRaw);
224     int batt = global.batteryState;
225
226     usnprintf(num, 30, "<— MEASURE DEBUG —>");
227     RIT128x96x4StringDraw(num, 0, 0, 15);
228
229     usnprintf(num, 30, "Raw temp: %d ", temp);
230     RIT128x96x4StringDraw(num, 0, 10, 15);
231
232     usnprintf(num, 30, "Raw Syst: %d ", sys);
233     RIT128x96x4StringDraw(num, 0, 20, 15);
234
235     usnprintf(num, 30, "Raw Dia: %d ", dia);
236     RIT128x96x4StringDraw(num, 0, 30, 15);
237
238     usnprintf(num, 30, "Raw Pulse: %d ", pulse);
239     RIT128x96x4StringDraw(num, 0, 40, 15);
240
241     usnprintf(num, 30, "Raw Batt: %d ", batt);
242     RIT128x96x4StringDraw(num, 0, 50, 15);

```

```

243 #endif
244
245     vTaskResume(computeHandle); // run the compute task
246 }

```

#### D.7.4 Compute Task

../code/compute.h

```

1  /*
2  * compute.h
3  * Author(s): PatrickMa
4  * 1/28/2014
5  *
6  * Defines the public interface for computeTask
7  * initializeComputeData() should be called before running computeTask()
8  */
9
10 #include "tcb.h"
11
12 /* Points to the TCB for compute */
13 extern TCB computeTask;

```

../code/compute.c

```

1  /*
2  * compute.c
3  * Author(s): PatrickMa
4  * 1/28/2014
5  *
6  * Implements compute.c
7  */
8
9  #define DEBUG_COMPUTE 0
10
11 #include "CircularBuffer.h"
12 #include "compute.h"
13 #include "globals.h"
14 #include "timebase.h"
15 #include "schedule.h"
16
17 // Used for debug display
18 #if DEBUG_COMPUTE
19 #include "drivers/rit128x96x4.h"
20 #include "utils/ustdlib.h"
21 #include "ekgProcess.h"
22 #endif
23
24 // computeData structure internal to compute task
25 typedef struct computeData {
26     // raw data pointers
27     CircularBuffer *temperatureRaw;

```

```

28 CircularBuffer *systolicPressRaw;
29 CircularBuffer *diastolicPressRaw;
30 CircularBuffer *pulseRateRaw;
31
32 //corrected data pointers
33 CircularBuffer *temperatureCorrected;
34 CircularBuffer *systolicPressCorrected;
35 CircularBuffer *diastolicPressCorrected;
36 CircularBuffer *pulseRateCorrected;
37
38 tBoolean *measurementComplete;
39 unsigned short *measurementSelect;
40 tBoolean *ekgCaptureDone;
41 tBoolean *ekgProcessDone;
42 } ComputeData;
43
44 void computeRunFunction(void *computeData);
45
46 static ComputeData data; // the internal data
47 TCB computeTask = {&computeRunFunction, &data}; // task interface
48
49 /*
50 * Initializes the computeData task values (pointers to variables, etc)
51 */
52 void initializeComputeTask() {
53     data.temperatureRaw = &(global.temperatureRaw);
54     data.systolicPressRaw = &(global.systolicPressRaw);
55     data.diastolicPressRaw = &(global.diastolicPressRaw);
56     data.pulseRateRaw = &(global.pulseRateRaw);
57
58     data.temperatureCorrected = &(global.temperatureCorrected);
59     data.systolicPressCorrected = &(global.systolicPressCorrected);
60     data.diastolicPressCorrected = &(global.diastolicPressCorrected);
61     data.pulseRateCorrected = &(global.pulseRateCorrected);
62
63     data.measurementComplete = &(global.measurementComplete);
64     data.measurementSelect = &(global.measurementSelection);
65     data.ekgCaptureDone = &(global.ekgCaptureDone);
66     data.ekgProcessDone = &(global.ekgProcessDone);
67 }
68
69 /*
70 * Linearizes the raw data measurement and converts value into human
71 * readable format
72 */
73 void computeRunFunction(void *computeData) {
74     static tBoolean onFirstRun = true;
75
76     if (onFirstRun) {
77         initializeComputeTask();
78         onFirstRun = false;
79     }
80

```



```

81 ComputeData *cData = (ComputeData *) computeData;
82
83 float temp = (5 + 0.75 * (*(int*)cbGet(cData->temperatureRaw)))/10;
84 float systolic = 9 + 2 * (*(int*)cbGet(cData->systolicPressRaw));
85 float diastolic = 6 + 1.5 * (*(int*)cbGet(cData->diastolicPressRaw));
86 float pulseRate = 8 + 3 * (*(int*)cbGet(cData->pulseRateRaw));
87
88 cbAdd(cData->temperatureCorrected, &temp);
89 cbAdd(cData->systolicPressCorrected, &systolic);
90 cbAdd(cData->diastolicPressCorrected, &diastolic);
91 cbAdd(cData->pulseRateCorrected, &pulseRate);
92
93 // if (0 == *(cData->measurementSelect) || 4 == *(cData->measurementSelect)
94 // {
95 //     while (!(cData->ekgCaptureDone)) { // wait until ekg captured
96 //     }
97 //     vTaskResume(ekgProcessHandle);
98 //     RIT128x96x4StringDraw("go ekgProcess", 0, 60, 15);
99 //     *(cData->ekgProcessDone) = false;
100 //     while (*(cData->ekgProcessDone)) { // wait until ekg computed
101 //     }
102 // }
103 *(cData->measurementComplete) = true;
104
105 vTaskSuspend(NULL); // suspend self
106
107 #if DEBUG_COMPUTE
108 char num[30];
109
110 usnprintf(num, 30, "<— COMPUTE DEBUG —>");
111 RIT128x96x4StringDraw(num, 0, 0, 15);
112
113 usnprintf(num, 30, "Corrected temp: %d", (unsigned int) temp);
114 RIT128x96x4StringDraw(num, 0, 10, 15);
115
116 usnprintf(num, 30, "Corrected Syst: %d", (unsigned int) systolic);
117 RIT128x96x4StringDraw(num, 0, 20, 15);
118
119 usnprintf(num, 30, "Corrected Dia: %d", (unsigned int) diastolic);
120 RIT128x96x4StringDraw(num, 0, 30, 15);
121
122 usnprintf(num, 30, "Corrected Pulse: %d", (unsigned int) pulseRate);
123 RIT128x96x4StringDraw(num, 0, 40, 15);
124 #endif
125 }

```

### D.7.5 Keypad Task

../code/keypad.h

```
1 /*
```

```

2  * keyPad.h
3  * Author(s): Jarrett Gaddy
4  * 2/10/2014
5  *
6  * Defines the public interface for the keyPad task
7  *
8  * initializeKeyPadTask() should be called once before performing runkeyPad
9  * functions
10 */
11
12 #include "tcb.h" // for TCBs
13
14 /* Initialize KeyPadData, must be done before running functions */
15 void initializeKeyPadTask();
16
17 /* The keyPad Task */
18 extern TCB keyPadTask;

```

../code/keypad.c

```

1  /*
2  * keypad.c
3
4  * Author(s): Jarrett Gaddy
5  * 2/10/2014
6  *
7  * implements keyPad.h
8  */
9
10 #include "keyPad.h"
11 #include "globals.h"
12 #include "timebase.h"
13 #include "inc/hw_types.h"
14 #include "driverlib/sysctl.h"
15 #include "driverlib/gpio.h"
16 #include "inc/hw_memmap.h"
17
18 #define UP_SW GPIOPinRead(GPIO_PORTD_BASE, GPIO_PIN_4)
19 #define DOWN_SW GPIOPinRead(GPIO_PORTD_BASE, GPIO_PIN_5)
20 #define LEFT_SW GPIOPinRead(GPIO_PORTD_BASE, GPIO_PIN_6)
21 #define RIGHT_SW GPIOPinRead(GPIO_PORTD_BASE, GPIO_PIN_7)
22 #define ACK_SW GPIOPinRead(GPIO_PORTE_BASE, GPIO_PIN_0)
23
24 // StatusData structure internal to compute task
25 typedef struct {
26     unsigned short *mode;
27     unsigned short *measurementSelection;
28     unsigned short *scroll;
29     tBoolean *alarmAcknowledge;
30     tBoolean *select;
31 } KeyPadData;
32
33 void keyPadRunFunction(void *data); // prototype for compiler

```

```

34
35 static KeyPadData data; // the internal data
36 TCB keyPadTask = {&keyPadRunFunction, &data}; // task interface
37
38 /* Initialize the StatusData task values */
39 void initializeKeyPadTask() {
40
41     SysCtlPeripheralEnable(SYSCTL_PERIPH_GPIOD);
42     SysCtlPeripheralEnable(SYSCTL_PERIPH_GPIOE);
43
44
45
46     GPIOPadConfigSet(GPIO_PORTE_BASE, GPIO_PIN_0, GPIO_STRENGTH_2MA,
47         GPIO_PIN_TYPE_STD_WPU);
48     GPIODirModeSet(GPIO_PORTE_BASE, GPIO_PIN_0, GPIO_DIR_MODE_IN);
49
50     GPIOPadConfigSet(GPIO_PORTD_BASE, GPIO_PIN_4, GPIO_STRENGTH_2MA,
51         GPIO_PIN_TYPE_STD);
52     GPIODirModeSet(GPIO_PORTD_BASE, GPIO_PIN_4, GPIO_DIR_MODE_IN);
53
54     GPIOPadConfigSet(GPIO_PORTD_BASE, GPIO_PIN_5, GPIO_STRENGTH_2MA,
55         GPIO_PIN_TYPE_STD);
56     GPIODirModeSet(GPIO_PORTD_BASE, GPIO_PIN_5, GPIO_DIR_MODE_IN);
57
58     GPIOPadConfigSet(GPIO_PORTD_BASE, GPIO_PIN_6, GPIO_STRENGTH_2MA,
59         GPIO_PIN_TYPE_STD);
60     GPIODirModeSet(GPIO_PORTD_BASE, GPIO_PIN_6, GPIO_DIR_MODE_IN);
61
62     GPIOPadConfigSet(GPIO_PORTD_BASE, GPIO_PIN_7, GPIO_STRENGTH_2MA,
63         GPIO_PIN_TYPE_STD);
64     GPIODirModeSet(GPIO_PORTD_BASE, GPIO_PIN_7, GPIO_DIR_MODE_IN);
65
66     //for ack switch
67     /* GPIOPadConfigSet(GPIO_PORTE_BASE, GPIO_PIN_3, GPIO_STRENGTH_2MA,
68         GPIO_PIN_TYPE_STD_WPU);
69     GPIODirModeSet(GPIO_PORTE_BASE, GPIO_PIN_3, GPIO_DIR_MODE_IN); */
70
71
72
73 // Load data
74 data.mode = &(global.mode);
75 data.measurementSelection = &(global.measurementSelection);
76 data.scroll = &(global.scroll);
77 data.alarmAcknowledge = &(global.alarmAcknowledge);
78 data.select = &(global.select);
79
80 // Load TCB
81 keyPadTask.runTaskFunction = &keyPadRunFunction;
82 keyPadTask.taskDataPtr = &data;
83 }
84
85 /* Perform status tasks */
86 void keyPadRunFunction(void *keyPadData){

```

```

87
88 static tBoolean onFirstRun = true;
89
90 if (onFirstRun) {
91     initializeKeyPadTask();
92     onFirstRun = false;
93 }
94 KeyPadData *kdata = (KeyPadData *) keyPadData;
95 if ( 0 == *(kdata->mode))
96 {
97
98     if (false == *(kdata->select))
99     {
100
101         if (UP_SW)
102             *(kdata->scroll) = *(kdata->scroll) + 1;
103         else if (DOWN_SW)
104             *(kdata->scroll) = *(kdata->scroll) - 1;
105         else if (RIGHT_SW)
106         {
107             *(kdata->select) = true;
108             *(kdata->measurementSelection) = (*(kdata->scroll) % 5) + 1;
109         }
110
111         else if (LEFT_SW)
112         {
113             *(kdata->mode) = 1;
114         }
115     }
116 }
117 else
118 {
119     if (LEFT_SW)
120     {
121         *(kdata->select) = false;
122         *(kdata->measurementSelection) = 0;
123     }
124 }
125 }
126 else
127 {
128     if (RIGHT_SW)
129     {
130         *(kdata->mode) = 0;
131     }
132 }
133 if (1 == ACK_SW)
134 {
135     *(kdata->alarmAcknowledge) = true;
136 }
137 else
138     *(kdata->alarmAcknowledge) = false;
139 }

```

### D.7.6 Display Task

../code/display.h

```
1 /*
2  * display.h
3  * Author(s): Jarrett Gaddy
4  * 1/28/2014, updated 2/10/2014
5  *
6  * Defines the interface for the displayTask.
7  * initializeDisplayData() should be called before running displayTask()
8  */
9
10 #include "tcb.h"
11
12 /* Initialize DisplayData, must be done before running displayTask() */
13 void initializeDisplayTask();
14
15 /* Points to the TCB for display */
16 extern TCB displayTask;
```

../code/display.c

```
1 /*
2  * display.c
3  * Author(s): jarrett Gaddy
4  * 2/10/2014
5  *
6  * Implements display.h
7  */
8
9 #include "globals.h"
10 #include "timebase.h"
11 #include "display.h"
12 #include "inc/hw_types.h"
13 #include "drivers/rit128x96x4.h"
14 #include "utils/ustdlib.h"
15 #include <stdlib.h>
16
17 #define DISPLAY_OFF 0
18
19 // Internal data structure
20 typedef struct oledDisplayData {
21     CircularBuffer *temperatureCorrected;
22     CircularBuffer *systolicPressCorrected;
23     CircularBuffer *diastolicPressCorrected;
24     CircularBuffer *pulseRateCorrected;
25     CircularBuffer *ekgFrequencyResult;
26     unsigned short *batteryState;
27     unsigned short *mode;
28     unsigned short *measurementSelection;
```

```

29 unsigned short *scroll;
30 tBoolean *alarmAcknowledge;
31 tBoolean *select;
32 unsigned short *selection;
33 } DisplayData;
34
35 void displayRunFunction(void *dataptr); // prototype for compiler
36
37 static DisplayData data; // internal data
38 TCB displayTask = {&displayRunFunction, &data}; // task interface
39
40 void initializeDisplayTask() {
41     // Load data
42     data.temperatureCorrected = &(global.temperatureCorrected);
43     data.systolicPressCorrected = &(global.systolicPressCorrected);
44     data.diastolicPressCorrected = &(global.diastolicPressCorrected);
45     data.pulseRateCorrected = &(global.pulseRateCorrected);
46     data.batteryState = &(global.batteryState);
47     data.ekgFrequencyResult = &(global.ekgFrequencyResult);
48
49     data.mode = &(global.mode);
50     data.measurementSelection = &(global.measurementSelection);
51     data.scroll = &(global.scroll);
52     data.alarmAcknowledge = &(global.alarmAcknowledge);
53     data.select = &(global.select);
54 }
55
56
57 void displayRunFunction(void *dataptr) {
58     static tBoolean onFirstRun = true;
59
60     if (onFirstRun) {
61         initializeDisplayTask();
62         onFirstRun = false;
63     }
64
65     DisplayData *dData = (DisplayData *) dataptr;
66
67     tBoolean selection = *(dData->select);
68     int scroll = *(dData->scroll);
69
70     #if !DISPLAY_OFF
71     char num[40];
72     //char buf1[30];
73     char buf2[30];
74
75     if (0 == *(dData->mode))
76     {
77         if (false == selection)
78         {
79             RIT128x96x4StringDraw("Make Selection", 0, 0,
15);

```

```

80     RIT128x96x4StringDraw( "   Blood Pressure           ", 0, 10,
15);
81     RIT128x96x4StringDraw( "   Temperature           ", 0, 20,
15);
82     RIT128x96x4StringDraw( "   Pulse Rate           ", 0, 30,
15);
83     RIT128x96x4StringDraw( "   EKG           ", 0, 40,
15);
84     RIT128x96x4StringDraw( "   Battery           ", 0, 50,
15);
85     RIT128x96x4StringDraw( "           ", 0, 60,
15);
86     RIT128x96x4StringDraw( "           ", 0, 70,
15);
87     RIT128x96x4StringDraw( "→", 0, 10*((scroll%5+1)), 15);
88 }
89 else
90 {
91     RIT128x96x4StringDraw( "           ", 0, 0,
15);
92     if (0 == scroll%5)
93         RIT128x96x4StringDraw( "Blood Pressure:", 0, 0, 15);
94     else if (1 == scroll%5)
95         RIT128x96x4StringDraw( "Temperature:", 0, 0, 15);
96     else if (2 == scroll%5)
97         RIT128x96x4StringDraw( "Pulse Rate:", 0, 0, 15);
98     else if (3 == scroll%5)
99         RIT128x96x4StringDraw( "EKG:", 0, 0, 15);
100    else if (4 == scroll%5)
101        RIT128x96x4StringDraw( "Battery:", 0, 0, 15);
102
103        else RIT128x96x4StringDraw( "oops", 0, 0, 15); //just in case
104
105
106
107        if (0 == scroll%5)
108            usnprintf(buf2,30, "Systolic: %d mm Hg ", (int) * ( (float*) cbGet(
dData->systolicPressCorrected)));
109        else if (1 == scroll%5)
110            usnprintf(buf2,30, "%d C ", (int) * ( (float*) cbGet(dData->
temperatureCorrected)));
111        else if (2 == scroll%5)
112            usnprintf(buf2,30, "%d BPM ", (int) * ( (float*) cbGet(dData->
pulseRateCorrected)));
113        else if (3 == scroll%5)
114            usnprintf(buf2,30, "%d Hz ", (int) * ( (float*) cbGet(dData->
ekgFrequencyResult)));
115        else if (4 == scroll%5)
116            usnprintf(buf2,30, "%d %% ", (int) *(dData->batteryState)/2);
117        //else buf2 = "oops"; //just in case
118
119        RIT128x96x4StringDraw( "           ", 0, 10,
15);

```

```

120     RIT128x96x4StringDraw(buf2, 0, 10, 15);
121     if(0 == scroll%5)
122     {
123         usnprintf(buf2,30, "Diastolic: %d mm Hg", (int)*(float*)
cbGet(dData->diastolicPressCorrected));
124
125         RIT128x96x4StringDraw(buf2, 0, 20, 15);
126     }
127     else RIT128x96x4StringDraw(" ",
0, 20, 15);
128     RIT128x96x4StringDraw(" ", 0, 30,
15);
129     RIT128x96x4StringDraw(" ", 0, 40,
15);
130     RIT128x96x4StringDraw(" ", 0, 50,
15);
131     RIT128x96x4StringDraw(" ", 0, 60,
15);
132     RIT128x96x4StringDraw(" ", 0, 70,
15);
133
134 }
135 }
136 else //(1 == *(data->mode))
137 {
138     usnprintf(num,40, "Temperature: %d C", (int) *(float*) cbGet
(dData->temperatureCorrected));
139     RIT128x96x4StringDraw(num, 0, 0, 15);
140
141     usnprintf(num,40, "Systolic Pressure: ");
142     RIT128x96x4StringDraw(num, 0, 10, 15);
143
144     usnprintf(num,40, "%d mm Hg", (int) *(float*) cbGet
(dData->systolicPressCorrected));
145     RIT128x96x4StringDraw(num, 0, 20, 15);
146
147     usnprintf(num,40, "Diastolic Pressure: ");
148     RIT128x96x4StringDraw(num, 0, 30, 15);
149
150     usnprintf(num,40, "%d mm Hg", (int) *(float*)
cbGet(dData->diastolicPressCorrected));
151     RIT128x96x4StringDraw(num, 0, 40, 15);
152
153     usnprintf(num,40, "Pulse rate: %d BPM", (int) *(float*)
cbGet(dData->pulseRateCorrected));
154     RIT128x96x4StringDraw(num, 0, 50, 15);
155
156     usnprintf(num,40, "EKG: %d Hz", *((int *) cbGet(
dData->ekgFrequencyResult));
157     RIT128x96x4StringDraw(num, 0, 60, 15);
158
159     usnprintf(num,40, "Battery: %d %%", (int) *(dData->
batteryState)/2);

```



```

160     RIT128x96x4StringDraw(num,0 , 70,15);
161 }
162 #endif
163 }

```

### D.7.7 Warning/Alarm Task

../code/warning.h

```

1  /*
2  * warning.h
3  * Author(s): Jarrett Gaddy
4  * 1/28/2014
5  *
6  * Defines the interface for the warning task
7  * initializeWarningData() should be called before running warningTask()
8  */
9
10 #include "tcb.h"
11
12 #define WARN_LOW 0.95 //warn at 5% below min range value
13 #define WARN_HIGH 1.05 //warn at 5% above max range value
14 #define ALARM_LOW 0.90 //alarm at 10% below min range value
15 #define ALARM_HIGH 1.20 //alarm at 20% above max range value
16
17
18 #define WARN_RATE_PULSE 4 // flash rate in terms of minor cycles
19 #define WARN_RATE_TEMP 2
20 #define WARN_RATE_PRESS 1
21
22 #define TEMP_MIN 36.1
23 #define TEMP_MAX 37.8
24 #define SYS_MAX 120
25 #define DIA_MAX 80
26 #define PULSE_MIN 60
27 #define PULSE_MAX 100
28 #define BATTERY_MIN 40
29
30 /* Initialize displayData, must be done before running warningTask() */
31 void initializeWarningTask();
32
33 /* The warning task */
34 extern TCB warningTask;

```

../code/warning.c

```

1  /*
2  * warning.c
3  * Author(s): jarrett Gaddy, PatrickMa
4  * 1/28/2014
5  *
6  * Implements warning.h

```

```

7  */
8
9  #define DEBUG_WARNING 0
10 #define ALARM_OFF 0
11
12 #include "globals.h"
13 #include "timebase.h"
14 #include "warning.h"
15 #include "schedule.h"
16 #include "inc/hw_types.h"
17 #include <stdlib.h>
18 #include <stdio.h>
19
20 #include "inc/hw_memmap.h"
21 #include "driverlib/debug.h"
22 #include "driverlib/gpio.h"
23 #include "driverlib/pwm.h"
24 #include "driverlib/sysctl.h"
25
26 #if DEBUG_WARNING
27 #include "drivers/rit128x96x4.h"
28 #include "utils/ustdlib.h"
29 #endif
30
31 // alarm cycle period (on/off) in millisecond
32 #define ALARM_PERIOD 2000
33 #define ALARM_CYCLE_RATE (ALARM_PERIOD / MINOR_CYCLE / 2)
34
35 // duration to sleep in terms of minor cycles
36 #define ALARM_SLEEP_PERIOD (5 * MAJOR_CYCLE)
37
38 #define LED_GREEN GPIO_PIN_6
39 #define LED_RED GPIO_PIN_5
40 #define LED_YELLOW GPIO_PIN_7
41
42 typedef enum {OFF, ON, ASLEEP} alarmState;
43 typedef enum {NONE, WARN_PRESS, WARN_TEMP, WARN_PULSE} warningState;
44 typedef enum {NORMAL, LOW} batteryState;
45
46 //pin E0 for input on switch 3
47 //pin C5 C6 and C7 for led out
48
49 // compiler prototypes
50 void warningRunFunction(void *dataptr); // prototype for compiler
51
52 // Internal data structure
53 typedef struct WarningData {
54     CircularBuffer *temperatureCorrected;
55     CircularBuffer *systolicPressCorrected;
56     CircularBuffer *diastolicPressCorrected;
57     CircularBuffer *pulseRateCorrected;
58     unsigned short *batteryState;
59 } WarningData;

```

```

60
61 extern tBoolean serialActive;
62 static unsigned long ulPeriod; // sets the alarm tone period
63 static WarningData data;      // internal data
64
65 TCB warningTask = {&warningRunFunction, &data}; // task interface
66
67
68 /*
69  * initializes task variables
70  */
71 void initializeWarningTask() {
72     //
73     // Enable the peripherals used by this code. I.e enable the use of pin
74     // banks, etc.
75     //
76     SysCtlPeripheralEnable(SYSCTL_PERIPH_PWM0);
77     SysCtlPeripheralEnable(SYSCTL_PERIPH_GPIOC); // bank C
78     SysCtlPeripheralEnable(SYSCTL_PERIPH_GPIOE); // bank E
79     SysCtlPeripheralEnable(SYSCTL_PERIPH_GPIOF); // bank F
80     SysCtlPeripheralEnable(SYSCTL_PERIPH_GPIOG); // bank G
81
82     // configure the pin C5 for 4mA output
83     GPIOPadConfigSet(GPIO_PORTC_BASE, LED_RED, GPIO_STRENGTH_4MA,
84                     GPIO_PIN_TYPE_STD);
85     GPIODirModeSet(GPIO_PORTC_BASE, LED_RED, GPIO_DIR_MODE_OUT);
86
87     // configure the pin C6 for 4mA output
88     GPIOPadConfigSet(GPIO_PORTC_BASE, LED_GREEN, GPIO_STRENGTH_4MA,
89                     GPIO_PIN_TYPE_STD);
90     GPIODirModeSet(GPIO_PORTC_BASE, LED_GREEN, GPIO_DIR_MODE_OUT);
91
92     // configure the pin C7 for 4mA output
93     GPIOPadConfigSet(GPIO_PORTC_BASE, LED_YELLOW, GPIO_STRENGTH_4MA,
94                     GPIO_PIN_TYPE_STD);
95     GPIODirModeSet(GPIO_PORTC_BASE, LED_YELLOW, GPIO_DIR_MODE_OUT);
96
97     // configure the pin E0 for input (sw3). NB: requires pull-up to operate
98     GPIOPadConfigSet(GPIO_PORTC_BASE, GPIO_PIN_0, GPIO_STRENGTH_2MA,
99                     GPIO_PIN_TYPE_STD_WPU);
100    GPIODirModeSet(GPIO_PORTC_BASE, GPIO_PIN_0, GPIO_DIR_MODE_IN);
101
102    /* This function call does the same result of the above pair of calls,
103     * but still requires that the bank of peripheral pins is enabled via
104     * SysCtlPeripheralEnable()
105     */
106    // GPIOPinTypeGPIOOutput(GPIO_PORTC_BASE, LED_RED);
107
108    // //////////////////////////////////////
109    // This section defines the PWM speaker characteristics
110    // //////////////////////////////////////

```

```

109
110 //
111 // Set the clocking to run directly from the crystal.
112 //
113 SysCtlPWMClockSet(SYSCTL_PWMDIV_1);
114
115 //
116 // Set GPIO F0 and G1 as PWM pins. They are used to output the PWM0 and
117 // PWM1 signals.
118 //
119 GPIOPinTypePWM(GPIO_PORTF_BASE, GPIO_PIN_0);
120 GPIOPinTypePWM(GPIO_PORTG_BASE, GPIO_PIN_1);
121
122 //
123 // Compute the PWM period based on the system clock.
124 //
125 uIPeriod = SysCtlClockGet() / 65;
126
127 //
128 // Set the PWM period to 440 (A) Hz.
129 //
130 PWMGenConfigure(PWM0_BASE, PWM_GEN_0,
131                 PWM_GEN_MODE_UP_DOWN | PWM_GEN_MODE_NO_SYNC);
132 PWMGenPeriodSet(PWM0_BASE, PWM_GEN_0, uIPeriod);
133
134 //
135 // Set PWM0 to a duty cycle of 25% and PWM1 to a duty cycle of 75%.
136 //
137 PWMPulseWidthSet(PWM0_BASE, PWM_OUT_0, uIPeriod / 4);
138 PWMPulseWidthSet(PWM0_BASE, PWM_OUT_1, uIPeriod * 3 / 4);
139
140 //
141 // Enable the PWM0 and PWM1 output signals.
142 //
143 PWMOutputState(PWM0_BASE, PWM_OUT_0_BIT | PWM_OUT_1_BIT, true);
144
145 // initialize the warning data pointers
146 data.temperatureCorrected = &(global.temperatureCorrected);
147 data.systolicPressCorrected = &(global.systolicPressCorrected);
148 data.diastolicPressCorrected = &(global.diastolicPressCorrected);
149 data.pulseRateCorrected = &(global.pulseRateCorrected);
150 data.batteryState = &(global.batteryState);
151
152 }
153
154 ///////////////////////////////////////////////////////////////////
155
156 /*
157  * Warning task function
158  */
159 void warningRunFunction(void *dataptr) {
160
161     static alarmState aState = OFF;

```

```

162 static warningState wState = NONE;
163 static batteryState bState = NORMAL;
164
165 static warningState prevState;
166 prevState = wState;
167
168 static int wakeUpAlarmAt = 0;
169
170 static tBoolean onFirstRun = true;
171
172 if (onFirstRun){
173     initializeWarningTask();
174     onFirstRun = false;
175 }
176
177
178 // Get measurement data
179 WarningData *wData = (WarningData *) dataptr;
180 float temp = *( (float*) cbGet(wData->temperatureCorrected));
181 float sysPress = *( (float*) cbGet(wData->systolicPressCorrected));
182 float diaPress = *( (float*) cbGet(wData->diastolicPressCorrected));
183 float pulse = *( (float*) cbGet(wData->pulseRateCorrected));
184 unsigned short battery = *(wData->batteryState);
185
186 // Alarm condition
187 if ( (sysPress > SYS_MAX*ALARM_HIGH) ) { //|| // Commented lines = lab2
188     // (temp < TEMP_MIN*ALARM_LOW || temp > (TEMP_MAX*ALARM_HIGH)) ||
189     // (diaPress > DIA_MAX*ALARM_HIGH) ||
190     // (pulse < PULSE_MIN*ALARM_LOW || pulse > PULSE_MAX*ALARM_HIGH) ) {
191
192     // Should only turn alarm ON if it was previously OFF. If it is
193     // ASLEEP, shouldn't do anything.
194     if (aState == OFF) aState = ON;
195 }
196 else
197     aState = OFF;
198
199 // Warning Condition
200 if ( sysPress > SYS_MAX*ALARM_HIGH || diaPress > DIA_MAX*ALARM_HIGH )
201     wState = WARN_PRESS;
202 else if ( temp < TEMP_MIN*WARN_LOW || temp > TEMP_MAX*WARN_HIGH )
203     wState = WARN_TEMP;
204 else if ( pulse < PULSE_MIN*WARN_LOW || pulse > PULSE_MAX*WARN_HIGH )
205     wState = WARN_PULSE;
206 else
207     wState = NONE;
208
209 // Battery Condition
210 if (battery < BATTERY_MIN)
211     bState = LOW;
212
213 // Handle speaker, based on alarm state
214 static tBoolean pwmEnable = false;

```

```

215  switch (aState) {
216      case ON:
217
218  #if !ALARM_OFF
219      if (0 == (minor_cycle_ctr % ALARM_CYCLE_RATE)) { // toggle between on/
220  off
221          if (pwmEnable)
222              PWMGenEnable(PWM0_BASE, PWM_GEN_0);
223          else
224              PWMGenDisable(PWM0_BASE, PWM_GEN_0);
225          pwmEnable = !pwmEnable;
226      }
227  #endif
228      break;
229  case ASLEEP:
230      PWMGenDisable(PWM0_BASE, PWM_GEN_0);
231      break;
232  default: // OFF
233      PWMGenDisable(PWM0_BASE, PWM_GEN_0);
234      break;
235  }
236
237  // Handle warning cases
238  static int toggletime;
239  switch (wState) {
240      case WARN_PRESS:
241          GPIOPinWrite(GPIO_PORTC_BASE, LED_GREEN, 0X00);
242
243          if (wState != prevState) {
244              GPIOPinWrite(GPIO_PORTC_BASE, LED_RED, 0XFF);
245              toggletime = WARN_RATE_PRESS;
246          }
247          else if (0 == minor_cycle_ctr%toggletime) {
248              if (GPIOPinRead(GPIO_PORTC_BASE, LED_RED) == 0)
249                  GPIOPinWrite(GPIO_PORTC_BASE, LED_RED, 0XFF);
250              else
251                  GPIOPinWrite(GPIO_PORTC_BASE, LED_RED, 0X00);
252          }
253
254          break;
255      case WARN_TEMP:
256          GPIOPinWrite(GPIO_PORTC_BASE, LED_GREEN, 0X00);
257
258          if (wState != prevState) {
259              GPIOPinWrite(GPIO_PORTC_BASE, LED_RED, 0XFF);
260              toggletime = WARN_RATE_TEMP;
261          }
262          else if (0 == minor_cycle_ctr%toggletime) {
263              if (GPIOPinRead(GPIO_PORTC_BASE, LED_RED) == 0)
264                  GPIOPinWrite(GPIO_PORTC_BASE, LED_RED, 0XFF);
265              else
266                  GPIOPinWrite(GPIO_PORTC_BASE, LED_RED, 0X00);

```

```

267     }
268     break;
269 case WARN_PULSE:
270     GPIOPinWrite(GPIO_PORTC_BASE, LED_GREEN, 0X00);
271
272     if (wState != prevState) {
273         GPIOPinWrite(GPIO_PORTC_BASE, LED_RED, 0XFF);
274         toggleTime = WARN_RATE_PULSE;
275     }
276     else if (0 == minor_cycle_ctr % toggleTime) {
277         if (GPIOPinRead(GPIO_PORTC_BASE, LED_RED) == 0)
278             GPIOPinWrite(GPIO_PORTC_BASE, LED_RED, 0XFF);
279         else
280             GPIOPinWrite(GPIO_PORTC_BASE, LED_RED, 0X00);
281     }
282     break;
283 default: // NORMAL
284     GPIOPinWrite(GPIO_PORTC_BASE, LED_GREEN, 0xFF);
285     GPIOPinWrite(GPIO_PORTC_BASE, LED_RED, 0x00);
286     break;
287 }
288
289 // activate the remote terminal task if in ANY warn or alarm state
290 if (NONE == wState && OFF == aState)
291     vTaskSuspend(serialHandle);
292 else
293     vTaskResume(serialHandle);
294
295 // battery state indicator
296 if (bState == LOW) {
297     GPIOPinWrite(GPIO_PORTC_BASE, LED_YELLOW, 0xFF);
298     GPIOPinWrite(GPIO_PORTC_BASE, LED_GREEN, 0x00);
299 }
300 else
301     GPIOPinWrite(GPIO_PORTC_BASE, LED_YELLOW, 0x00);
302
303 /* This is the alarm override
304  * Upon override, the alarm is silenced for some time.
305  * silence length is defined by ALARM_SLEEP_PERIOD
306  *
307  * If the button is pushed, the value returned is 0
308  * If the button is NOT pushed, the value is non-zero
309  */
310 if ( ( 0 == global.alarmAcknowledge) && (aState == ON) )
311 {
312     // GPIOPinWrite(GPIO_PORTC_BASE, LED_YELLOW, 0XFF); // for debug, lights
313     // led
314     aState = ASLEEP;
315     wakeUpAlarmAt = minor_cycle_ctr + ALARM_SLEEP_PERIOD;
316 }
317
318 // Check whether to resound alarm
319 if (minor_cycle_ctr == wakeUpAlarmAt && aState == ASLEEP) {

```

```

319     aState = ON;
320     // GPIOPinWrite(GPIO_PORTC_BASE, LED_YELLOW, 0X00); // for debug, kills
    led
321 }
322
323 #if DEBUG_WARNING
324     char num[30];
325
326     usnprintf(num, 30, "<— WARNING DEBUG —>");
327     RIT128x96x4StringDraw(num, 0, 0, 15);
328
329     usnprintf(num, 30, "Cor temp: %d ", (int) temp);
330     RIT128x96x4StringDraw(num, 0, 10, 15);
331
332     usnprintf(num, 30, "Cor Syst: %d ", (int) sysPress);
333     RIT128x96x4StringDraw(num, 0, 20, 15);
334
335     usnprintf(num, 30, "Cor Dia: %d ", (int) diaPress);
336     RIT128x96x4StringDraw(num, 0, 30, 15);
337
338     usnprintf(num, 30, "Cor Pulse: %d ", (int) pulse);
339     RIT128x96x4StringDraw(num, 0, 40, 15);
340
341     usnprintf(num, 30, "Cor Batt: %d ", (unsigned short) battery);
342     RIT128x96x4StringDraw(num, 0, 50, 15);
343
344
345     usnprintf(num, 30, "aState: %d ", aState);
346     RIT128x96x4StringDraw(num, 0, 60, 15);
347
348     usnprintf(num, 30, "alarmAck: %d ", global.alarmAcknowledge);
349     RIT128x96x4StringDraw(num, 0, 70, 15);
350
351     usnprintf(num, 30, "pwmEn: %d ", pwmEnable);
352     RIT128x96x4StringDraw(num, 0, 80, 15);
353 #endif
354
355 }

```

### D.7.8 Serial Task

../code/serial.h

```

1  /*
2  *  serial.h
3  *  Author(s): Jonathan Ellington
4  *  2/05/2014
5  *
6  *  Defines the serial communications task. This sends various data over
7  *  RS-232.
8  *
9  *  This function should only run if there is a warning (in other words, the

```



```

10  * scheduler should add it to the task queue in the event of a warning), and
11  * should subsequently delete itself from the task queue.
12  */
13
14  #include "tcb.h" // for TCBs
15
16  /* The status Task */
17  extern TCB serialTask;

```

../code/serial.c

```

1  /*
2  * serial.c
3  * Author(s): Jonathan Ellington
4  * 2/05/2014
5  *
6  * Implements serial.c
7  */
8
9  #include "globals.h"
10 #include "timebase.h"
11 #include "serial.h"
12 #include "schedule.h"
13 #include "CircularBuffer.h"
14 #include "inc/hw_types.h"
15 #include "driverlib/uart.h"
16 #include "driverlib/gpio.h"
17 #include "inc/hw_memmap.h"
18 #include "utils/ustdlib.h"
19 #include <string.h>
20
21
22 // Internal data structure
23 typedef struct serialData {
24     CircularBuffer *temperatureCorrected;
25     CircularBuffer *systolicPressCorrected;
26     CircularBuffer *diastolicPressCorrected;
27     CircularBuffer *pulseRateCorrected;
28     CircularBuffer *ekgFrequencyResult;
29     unsigned short *batteryState;
30 } SerialData;
31
32 // Prototype
33 void UARTSend(const unsigned char *pucBuffer, unsigned long ulCount);
34 void serialRunFunction(void *dataptr);
35
36 static SerialData data; // internal data
37 TCB serialTask = {&serialRunFunction, &data}; // task interface
38
39 void initializeSerialTask() {
40     // Initialize Data
41     data.temperatureCorrected = &(global.temperatureCorrected);
42     data.systolicPressCorrected = &(global.systolicPressCorrected);

```

```

43 data.diastolicPressCorrected = &(global.diastolicPressCorrected);
44 data.pulseRateCorrected = &(global.pulseRateCorrected);
45 data.batteryState = &(global.batteryState);
46 data.ekgFrequencyResult = &(global.ekgFrequencyResult);
47 // UART Stuff
48 // Enable the peripherals used by this example.
49 SysCtlPeripheralEnable(SYSCTL_PERIPH_UART0);
50 SysCtlPeripheralEnable(SYSCTL_PERIPH_GPIOA);
51
52 // Set GPIO A0 (UART RX)
53 GPIOPinTypeUART(GPIO_PORTA_BASE, GPIO_PIN_1);
54
55 // Configure the UART for 115,200, 8–N–1 operation.
56 UARTConfigSetExpClk(UART0_BASE, SysCtlClockGet(), 115200,
57     (UART_CONFIG_WLEN_8 | UART_CONFIG_STOP_ONE |
58     UART_CONFIG_PAR_NONE));
59
60 UARTEnable(UART0_BASE);
61 }
62
63 void serialRunFunction(void *dataptr) {
64     static tBoolean onFirstRun = true;
65     SerialData *sData = (SerialData *) dataptr;
66     if (onFirstRun) {
67         initializeSerialTask();
68         onFirstRun = false;
69     }
70
71
72     int temp = (int) *((float *)cbGet(sData->temperatureCorrected));
73     int sys = (int) *((float *)cbGet(sData->systolicPressCorrected));
74     int dia = (int) *((float *)cbGet(sData->diastolicPressCorrected));
75     int pulse = (int) *((float *)cbGet(sData->pulseRateCorrected));
76     int EKG = *((int *)cbGet(sData->ekgFrequencyResult));
77     int batt = *(data.batteryState);
78
79     char buf[40];
80     usnprintf(buf, 40, "\f1. Temperature:\t\t%d C\n\nr", temp);
81     UARTSend( (unsigned char *) buf, strlen(buf));
82
83     usnprintf(buf, 40, "2. Systolic pressure:\t\t%d mm Hg\n\nr", sys);
84     UARTSend( (unsigned char *) buf, strlen(buf));
85
86     usnprintf(buf, 40, "3. Diastolic pressure:\t\t%d mm Hg\n\nr", dia);
87     UARTSend( (unsigned char *) buf, strlen(buf));
88
89     usnprintf(buf, 40, "4. Pulse rate:\t\t%d BPM\n\nr", pulse);
90     UARTSend( (unsigned char *) buf, strlen(buf));
91
92     usnprintf(buf, 40, "5. EKG:\t\t\t%d Hz\n\nr", EKG);
93     UARTSend( (unsigned char *) buf, strlen(buf));
94
95     usnprintf(buf, 40, "6. Battery:\t\t\t%d%%\n\nr", batt/2);

```

```

96 UARTSend( (unsigned char *) buf, strlen(buf));
97
98 vTaskSuspend(NULL);    // suspend self
99 }
100
101 void UARTSend(const unsigned char *pucBuffer, unsigned long ulCount) {
102     // Loop while there are more characters to send.
103     while(ulCount--)
104     {
105         // Write the next character to the UART.
106         // This blocks while the FIFO queue is full
107         UARTCharPut(UART0_BASE, *pucBuffer++);
108     }
109 }

```

### D.7.9 EKG Capture

../code/ekgCapture.c

```

1  /* Author: patrick Ma
2   * 2/21/14
3   *
4   * ekgcapture.c
5   *
6   * Reads and stores data from ekg sensor via ADC. Stores data into a memory
7   * buffer.
8   */
9
10 #define DEBUG_EKG 0 // ekg task debug
11
12 #include "FreeRTOS.h"
13 #include "task.h"
14
15 #include "inc/hw_types.h"
16 #include "inc/hw_memmap.h"
17 #include "driverlib/adc.h" // for ADC use
18 #include "driverlib/timer.h" // for hw timer use
19 #include "driverlib/interrupt.h"
20 #include "inc/hw_ints.h"
21
22 #include "schedule.h"
23 #include "globals.h"
24 #include "ekgCapture.h"
25 #include "timebase.h"
26
27 // Used for debug display
28 #if DEBUG_EKG
29 #include "drivers/rit128x96x4.h"
30 #include "utils/ustdlib.h"
31 #endif
32
33 #if DEBUG_EKG

```

```

34 static char num[30]; // used for display
35 #endif
36
37 static tBoolean firstRun = true;
38 static tBoolean ekgComplete = false;
39 extern xTaskHandle ekgProcessHandle;
40
41 // ekgCapture data structure. Internal to the task
42 typedef struct ekgCaptureData {
43     signed int *ekgRawDataAddr; // raw output array address
44     tBoolean *ekgCaptureDone;
45 }EKGCaptureData;
46
47 static EKGCaptureData data; // the internal data
48
49 void ekgCaptureRunFunction(void *ekgCaptureData); // Compiler function
         prototypes
50 TCB ekgCaptureTask = {&ekgCaptureRunFunction, &data}; // task interface
51 void ADC0IntHandler();
52
53 /*
54  * sets up task specific variables, etc
55  */
56 void initializeEKGTask() {
57     #if DEBUG_EKG
58         RIT128x96x4Init(1000000);
59         RIT128x96x4StringDraw("EKG Capture Debug", 0, 0, 15);
60     #endif
61
62     data.ekgRawDataAddr = global.ekgRaw;
63     data.ekgCaptureDone = &(global.ekgCaptureDone);
64
65     // Enable read from GPIO pin (move this and below to startup?)
66     SysCtlPeripheralEnable(SYSCTL_PERIPH_GPIOE);
67     GPIOPinTypeADC(GPIO_PORTE_BASE, GPIO_PIN_7);
68
69     // Setup ADC0 for EKG capture
70     SysCtlPeripheralEnable(SYSCTL_PERIPH_ADC0);
71     SysCtlADCSpeedSet(SYSCTL_ADCSPEED_500KSPS);
72     ADCSequenceDisable(ADC0_BASE, EKG_SEQ);
73     ADCSequenceConfigure(ADC0_BASE, EKG_SEQ, ADC_TRIGGER_PROCESSOR, 1);
74     ADCSequenceStepConfigure(ADC0_BASE, EKG_SEQ, 0, ADC_CTL_IE | ADC_CTL_END |
        EKG.CH);
75     ADCIntRegister(ADC0_BASE, EKG_SEQ, ADC0IntHandler);
76     ADCIntEnable(ADC0_BASE, EKG_SEQ);
77
78     // configure timer (uses both 16-bit timers) for periodic timing
79     // SysCtlPeripheralEnable(EKG_TIMER_PERIPH);
80     // TimerDisable(EKG_TIMER_BASE, EKG_TIMER);
81     // TimerConfigure(EKG_TIMER_BASE, EKG_CFG);
82     // TimerControlTrigger(EKG_TIMER_BASE, EKG_TIMER, true);
83     // TimerLoadSet(EKG_TIMER_BASE, EKG_TIMER, (SysCtlClockGet() / SAMPLE_FREQ)
        );

```

```

84 // TimerEnable(EKG_TIMER_BASE, EKG_TIMER);
85
86 #if DEBUG_EKG
87     long clk = SysCtlClockGet();
88     long timeLoad = TimerLoadGet(EKG_TIMER_BASE, EKG_TIMER);
89     usnprintf(num, 30, "load: %ul", timeLoad);
90     RIT128x96x4StringDraw(num, 0, 10, 15);
91     usnprintf(num, 30, "clock: %ul", clk);
92     RIT128x96x4StringDraw(num, 0, 20, 15);
93 #endif
94 }
95
96 void delay_in_ms(int ms) {
97     for (volatile int i = 0; i < ms; i++)
98         for (volatile int j = 0; j < 500; j++);
99 }
100
101 /*
102  * captures a sequence of samples (given by NUM_EKG_SAMPLES) via the ADC and
103  * stores the results into a buffer
104  */
105 void ekgCaptureRunFunction(void *ekgCaptureData) {
106     if (firstRun) {
107         firstRun = false;
108         initializeEKGTask();
109     }
110
111     EKGCaptureData *eData = (EKGCaptureData *) ekgCaptureData;
112     ekgComplete = false; // reset the adc counters
113
114     ADCSequenceEnable(ADC0_BASE, EKG_SEQ);
115     while (!ekgComplete) {
116         ADCProcessorTrigger(ADC0_BASE, EKG_SEQ);
117         delay_in_ms(1);
118     }
119     ADCSequenceDisable(ADC0_BASE, EKG_SEQ);
120
121     *(eData->ekgCaptureDone) = true; // we want to process our measurement
122
123 #if DEBUG_EKG
124     usnprintf(num, 30, "end ADC get");
125     RIT128x96x4StringDraw(num, 0, 40, 15);
126
127     // let's check our values
128
129     int a = global.ekgRaw[0];
130     int b = global.ekgRaw[50];
131     int c = (int) *(data.ekgRawDataAddr + 255);
132     usnprintf(num, 30, "%d, %d, %d", a, b, c);
133     RIT128x96x4StringDraw(num, 3, 50, 15);
134 #endif
135
136     vTaskResume(ekgProcessHandle);

```

```

137 }
138
139 void ADC0IntHandler() {
140     unsigned long value;
141     static int i = 0;
142     //debugPin47();
143     // Read the value from the ADC.
144
145     while (1 != ADCSequenceDataGet(ADC0_BASE, EKG_SEQ, &value));
146     (data.ekgRawDataAddr)[i++] = (signed int) value;
147
148     // Done sampling if we've taken enough samples
149     if (i >= NUM.EKG_SAMPLES) {
150         i = 0;
151         ekgComplete = true;
152     }
153     ADCIntClear(ADC0_BASE, EKG_SEQ);
154 }

```

../code/ekgCapture.h

```

1 /*
2  * EkgCapture.h
3  * Author: Patrick Ma
4  * 2/22/2014
5  *
6  * Defines the public interface for the EKG data capture task.
7  *
8  */
9
10 #include "tcb.h"
11
12 #define EKG_TIMER_BASE TIMER0_BASE // the timer base used for sample
    collection
13 #define EKG_TIMER TIMER_B // the timer used for EKG sample
    collection
14 #define EKG_CFG TIMER_CFG_B_PERIODIC
15 #define EKG_TIMER_PERIPH SYSCTL_PERIPH_TIMER0
16
17 #define EKG_SEQ 0 // The sequence number assigned to the ekg sensor
18 #define EKG_CH ADC_CTL_CH0 // EKG analog input channel (ch0: pinE7, others:
    pinE4-6?)
19 #define EKG_PRIORITY 1 // the ekg sequence capture priority
20
21 /* Points to the TCB for the task */
22 extern TCB ekgCaptureTask;

```

#### D.7.10 EKG Process

../code/ekgProcess.c

```

1 /*

```

```

2  * ekgProcess.c
3  * Author: Patrick Ma
4  * 2/23/2014
5  *
6  * Processes the raw ekg samples and produces a frequency value
7  *
8  */
9
10 #define DEBUG_PROC 0
11
12 #include "FreeRTOS.h"
13 #include "task.h"
14
15 #include "globals.h"
16 #include "ekgProcess.h"
17 #include "CircularBuffer.h"
18 #include "optfft.h"
19 #include <stdio.h>
20 #include <string.h>
21
22 // Used for debug display
23 #if DEBUG_PROC
24 #include "drivers/rit128x96x4.h"
25 #include "utils/ustdlib.h"
26 static char num[30];
27 #endif
28
29 static tBoolean firstRun = true;
30
31 // ekgProcess data structure. Internal to this task
32 typedef struct egkProcessData {
33     signed int *ekgRawData;
34     signed int *ekgImgData;
35     CircularBuffer *ekgFreqResult;
36 } EKGProcessData;
37
38 static EKGProcessData data; // internal data object
39
40 void ekgProcessRunFunction(void *ekgProcessData);
41
42 TCB ekgProcessTask = {&ekgProcessRunFunction, &data}; // task interface
43
44 void initializeEKGProcess() {
45     data.ekgRawData = (global.ekgRaw);
46     data.ekgImgData = (global.ekgTemp);
47     data.ekgFreqResult = &(global.ekgFrequencyResult);
48
49 #if DEBUG_PROC
50     RIT128x96x4Init(1000000);
51     RIT128x96x4StringDraw("EKGProcess Debug", 0, 0, 15);
52 #endif
53 }
54

```

```

55 |
56 | /*
57 |  * Reads in the EKG samples and performs a fast Fourier transform (FFT) on
58 |  * the
59 |  * data to extract the primary frequency of the signal
60 |  */
61 | void ekgProcessRunFunction(void *ekgData) {
62 |     EKGProcessData *eData = (EKGProcessData *) ekgData;
63 |     if (firstRun) {
64 |         firstRun = false;
65 |         initializeEKGProcess();
66 |     }
67 |
68 |     // need to bit shift >> 4 (divide 16) then subtract 32
69 |     int i = 0;
70 |     int t = 0; // debug
71 |     for (i = 0; i < NUM_EKG_SAMPLES; i++) {
72 | #if DEBUG_PROC
73 |         usnprintf(num, 30, "%d \n", eData->ekgRawData[i]);
74 |         RIT128x96x4StringDraw(num, 0, 10, 15);
75 | #endif
76 |         int d = ((eData->ekgRawData[i] >> 4) - 31;
77 |         eData->ekgRawData[i] = d;
78 | #if DEBUG_PROC
79 |         if (eData->ekgRawData[i] > eData->ekgRawData[t])
80 |             t = i;
81 |         usnprintf(num, 30, "%d : %d \n", eData->ekgRawData[i], eData->ekgRawData
82 |             [t]);
83 |         RIT128x96x4StringDraw(num, 0, 20, 15);
84 | #endif
85 |     }
86 |
87 |     // reset Imaginary array
88 |     memset(eData->ekgImgData, 0, sizeof(signed int) * NUM_EKG_SAMPLES);
89 |     // memset(eData->ekgRawData, 0, sizeof(signed int) * NUM_EKG_SAMPLES);
90 |
91 |     signed int max_index = optfft( eData->ekgRawData, eData->ekgImgData );
92 |     //post processing
93 |     int freq = ((9166) * max_index) / NUM_EKG_SAMPLES;
94 |
95 | #if DEBUG_PROC
96 |     usnprintf(num, 30, "%d : %d ", max_index, freq);
97 |     RIT128x96x4StringDraw(num, 0, 30, 15);
98 | #endif
99 |     cbAdd(eData->ekgFreqResult, &freq);
100 |     vTaskSuspend(NULL);
101 | }

```

../code/ekgProcess.h

```

1 | /*
2 |  * ekgProcess.h

```



```

3  * Author: Patrick Ma
4  * 2/23/2014
5  *
6  * Header and public interface of the ekgProcess task
7  *
8  */
9
10 #include "tcb.h"
11
12 /* Points to the TCB for the task */
13
14 extern TCB ekgProcessTask;

```

### D.7.11 Command Task

../code/commandTask.c

```

1  /*
2  * commandTask.c
3  * Author(s); Patrick Ma
4  *
5  * 3/10/2014
6  *
7  * Interprets the text commands from the remote connection or system,
8  * performs
9  * the actions requested, and sends a reply signal.
10 */
11 #define DEBUG_COMMAND 0
12
13 #include "FreeRTOS.h"
14 #include "task.h"
15 #include "globals.h"
16 #include "warning.h"
17 #include "CircularBuffer.h"
18 #include "commandTask.h"
19 #include <string.h>
20 #include "hw_ints.h"
21 #include "driverlib/interrupt.h"
22 #include "driverlib/adc.h"
23
24
25 #include "drivers/rit128x96x4.h"
26 #include "utils/ustdlib.h"
27 char num[30];
28
29
30 #define TOKEN_DELIM " \t" // token delimiter values
31 #define TEMP_BUFFER_LEN 55
32
33 // compiler prototypes
34 void commandRunFunction(void *commandDataPtr);

```

```

35
36 // internal command data structure
37 typedef struct commandData
38 {
39     char *commandStr;
40     char *responseStr;
41     unsigned short *measureSelect;
42     tBoolean *measureComplete;
43     tBoolean *responseReady;
44     CircularBuffer *temperature;
45     CircularBuffer *systPress;
46     CircularBuffer *diasPress;
47     CircularBuffer *pulse;
48     CircularBuffer *ekg;
49     unsigned short *battery;
50 } CommandData;
51
52 static CommandData data; // version of data exposed to outside
53 TCB commandTask = {&commandRunFunction, &data}; // set up task interface
54 extern xTaskHandle measureHandle;
55 extern xTaskHandle displayHandle;
56 extern xTaskHandle ekgCaptureHandle;
57
58 /*
59  * local private variables
60  */
61 static tBoolean initialized = false;
62 static tBoolean measureOn;
63 static tBoolean displayOn = true;
64 static int value;
65 static char *cmd;
66 static char *sensor;
67
68 // communication arrays
69 static char parseArr[COMMAND_LENGTH];
70 static char temporaryBuffer[TEMP_BUFFER_LEN];
71 static char roughString[TEMP_BUFFER_LEN];
72 static char *formattedStr = temporaryBuffer;
73
74
75 /*
76  * initializes task variables
77  */
78 void initializeCommandTask(){
79     data.commandStr = (global.commandStr);
80     data.responseStr = (global.responseStr);
81     data.measureSelect = &(global.measurementSelection);
82     data.measureComplete = &(global.measurementComplete);
83     data.responseReady = &(global.responseReady);
84
85     data.temperature = &(global.temperatureCorrected);
86     data.systPress = &(global.systolicPressCorrected);
87     data.diasPress = &(global.diastolicPressCorrected);

```

```

88 data.pulse = &(global.pulseRateCorrected);
89 data.ekg = &(global.ekgFrequencyResult);
90 data.battery = &(global.batteryState);
91 }
92
93 /*
94  * parses the command string for command parameters, writing those values
95  * to
96  * variables
97  */
98 void parse(CommandData *cData) {
99     char delim[2] = TOKEN_DELIM;
100
101     strncpy(parseArr, cData->commandStr, COMMAND_LENGTH - 1);
102
103     cmd = strtok(parseArr, delim);
104     sensor = strtok(NULL, delim);
105 }
106
107 /*
108  * Adds an acknowledge or not acknowledge to the response buffer
109  */
110 void ackNack(CommandData *cData, tBoolean stat) {
111     if (stat)
112         strcat(cData->responseStr, "<p>A</p>", RESPONSE_LENGTH - 1);
113     else
114         strcat(cData->responseStr, "<p>E</p>", RESPONSE_LENGTH - 1);
115 }
116
117 /*
118  * Formats given string with appropriate html tags. Returns pointer to
119  * formatted string. If statusOK is false, <blink> or </blink> flags are
120  * added.
121  * CAUTION! addTags does not guard against buffer overflow. Make sure your
122  * string is not too long. The tags add up to 22 characters.
123  */
124 char* addTags(char* string, tBoolean statusOK) {
125     memset(temporaryBuffer, '\0', TEMP_BUFFER_LEN - 1);
126     strncpy(temporaryBuffer, "<p>", 3); // first tag
127
128     if (!statusOK) {
129         strcat(temporaryBuffer, "<blink>"); // include warning maybe
130     }
131     strcat(temporaryBuffer, string);
132     if (!statusOK) {
133         strcat(temporaryBuffer, "</blink>");
134     }
135     strcat(temporaryBuffer, "</p>");
136
137     return temporaryBuffer;
138 }

```

```

139 /*
140  * measures the data from a sensor
141  */
142 void measureFromSensor(CommandData* cData) {
143     tBoolean meas = true;
144     switch (*sensor) {
145         case 'A' :
146 #if DEBUG.COMMAND
147         RIT128x96x4StringDraw("take A", 2, 30, 15);
148 #endif
149         *(cData->measureSelect) = 0;
150         break;
151         case 'T' :
152 #if DEBUG.COMMAND
153         RIT128x96x4StringDraw("take T", 2, 30, 15);
154 #endif
155         *(cData->measureSelect) = 1;
156         break;
157         case 'B' :
158 #if DEBUG.COMMAND
159         RIT128x96x4StringDraw("take B", 2, 30, 15);
160 #endif
161         *(cData->measureSelect) = 2;
162         break;
163         case 'P' :
164 #if DEBUG.COMMAND
165         RIT128x96x4StringDraw("take P", 2, 30, 15);
166 #endif
167         *(cData->measureSelect) = 3;
168         break;
169         case 'E' :
170 #if DEBUG.COMMAND
171         RIT128x96x4StringDraw("take E", 2, 30, 15);
172 #endif
173         *(cData->measureSelect) = 4;
174         break;
175         case 'S' :
176 #if DEBUG.COMMAND
177         RIT128x96x4StringDraw("unsupported", 0, 30, 15);
178 #endif
179         *(cData->measureSelect) = '%';
180         break;
181         default :
182 #if DEBUG.COMMAND
183         RIT128x96x4StringDraw("invalid command", 0, 30, 15);
184 #endif
185         meas = false;
186     }
187
188     if (meas) {
189         ackNack(cData, true);
190 #if DEBUG.COMMAND
191         RIT128x96x4StringDraw("MeasureTask go!", 0, 40, 15);

```

```

192 #endif
193 } else {
194 #if DEBUG_COMMAND
195     RIT128x96x4StringDraw("no Measure", 0, 40, 15);
196 #endif
197     ackNack(cData, false);
198 }
199
200 // while (!(cData->measureComplete)) { // wait until measurement is
201 //     finished
202 // }
203
204 /*
205  * gets & formats string
206  */
207 void printTemp(CommandData *cData, tBoolean statusOK) {
208     value = (int) *(float *)cbGet(cData->temperature);
209     usnprintf(roughString, TEMP_BUFFER_LEN, "Temperature: %d C", value);
210     addTags(roughString, statusOK);
211     strncat(cData->responseStr, formattedStr, RESPONSE_LENGTH - strlen(cData->
212         responseStr) - 1);
213 }
214 /*
215  * gets & formats string
216  */
217 void printPressure(CommandData *cData, tBoolean statusOK) {
218     value = (int) *(float *)cbGet(cData->systPress);
219     usnprintf(roughString, TEMP_BUFFER_LEN, "Systolic Pressure: %d mmHg",
220         value);
221     addTags(roughString, statusOK);
222     strncat(cData->responseStr, formattedStr, RESPONSE_LENGTH - strlen(cData->
223         responseStr) - 1);
224     value = (int) *(float *)cbGet(cData->diasPress);
225     usnprintf(roughString, TEMP_BUFFER_LEN, "Diastolic Pressure: %d mmHg",
226         value);
227     addTags(roughString, statusOK);
228     strncat(cData->responseStr, formattedStr, RESPONSE_LENGTH - strlen(cData->
229         responseStr) - 1);
230 }
231 /*
232  * gets & formats string
233  */
234 void printPulse(CommandData *cData, tBoolean statusOK) {
235     value = (int) *(float *)cbGet(cData->pulse);
236     usnprintf(roughString, TEMP_BUFFER_LEN, "Pulse Rate: %d BPM", value);
237     addTags(roughString, statusOK);
238     strncat(cData->responseStr, formattedStr, RESPONSE_LENGTH - strlen(cData->
239         responseStr) - 1);
240 }

```

```

238 /*
239  * gets & formats string
240  */
241 void printEKG(CommandData *cData, tBoolean statusOK) {
242     value = *(int *)cbGet(cData->ekg);
243     usnprintf(roughString, TEMP_BUFFER_LEN, "EKG Frequency: %d Hz", value);
244     addTags(roughString, statusOK);
245     strncat(cData->responseStr, formattedStr, RESPONSE_LENGTH - strlen(cData->
        responseStr) - 1);
246 }
247
248 /*
249  * gets & formats string
250  */
251 void printBattery(CommandData *cData, tBoolean statusOK) {
252     value = (int) *(cData->battery) / 2;
253     usnprintf(roughString, TEMP_BUFFER_LEN, "Battery Remaining: %d %%", value)
        ;
254     addTags(roughString, statusOK);
255     strncat(cData->responseStr, formattedStr, RESPONSE_LENGTH - strlen(cData->
        responseStr) - 1);
256 }
257
258 /*
259  * Format the outgoing response string based on the specified measurement
260  */
261 void formatResponseStr(CommandData *cData){
262
263     int temp = *(int *)cbGet(cData->temperature);
264     int sysPress = *(int *)cbGet(cData->systPress);
265     int diaPress = *(int *)cbGet(cData->diasPress);
266     int pulse = *(int *)cbGet(cData->pulse);
267
268     tBoolean bpWarn = false;
269     tBoolean tempWarn = false;
270     tBoolean pulseWarn = false;
271     tBoolean battWarn = false;
272
273     if ( sysPress > SYS_MAX*ALARM_HIGH || diaPress > DIA_MAX*ALARM_HIGH )
274         bpWarn = true;
275     if ( temp < TEMP_MIN*WARN_LOW || temp > TEMP_MAX*WARN_HIGH )
276         tempWarn = true;
277     if ( pulse < PULSE_MIN*WARN_LOW || pulse > PULSE_MAX*WARN_HIGH )
278         pulseWarn = true;
279     if (*(cData->battery) < BATTERY_MIN)
280         battWarn = true;
281
282     switch (*sensor) {
283     case 'A' : //take all measurements
284         printTemp(cData, !tempWarn);
285         printPressure(cData, !bpWarn);
286         printPulse(cData, !pulseWarn);
287         printEKG(cData, true);

```

```

288     printBattery(cData, !battWarn);
289     break;
290 case 'T' : // get temperature measurementk
291     printTemp(cData, tempWarn);
292     break;
293 case 'B' : // get syst
294     printPressure(cData, bpWarn);
295     break;
296 case 'P' : //pulse rate
297     printPulse(cData, pulseWarn);
298     break;
299 case 'E' : //ekg frequency
300     printEKG(cData, true);
301     break;
302 case 'S' : // battery state
303     printBattery(cData, battWarn);
304     break;
305 }
306 }
307
308 /*
309  * runs the command task
310  *
311  * NB: the response string needs html formatting:
312  * ex: <p> Temperature: 50 C </p> <p> <blink> Blood Pressure: 120 </blink>
313  * </p>
314  */
314 void commandRunFunction(void *commandDataPtr) {
315     CommandData *cData = (CommandData *) commandDataPtr;
316     *(cData->responseReady) = false;
317     if (!initialized) {
318         initialized = true;
319         initializeCommandTask();
320     }
321
322 #if DEBUG.COMMAND
323     usnprintf(num, 30, "Initialize cmd function");
324     RIT128x96x4StringDraw(num, 0, 10, 15);
325 #endif
326
327     parse(cData);
328
329 #if DEBUG.COMMAND
330     usnprintf(num, 30, "%s %s", cmd, sensor);
331     RIT128x96x4StringDraw(num, 0, 20, 15);
332 #endif
333
334     memset(cData->responseStr, '\0', RESPONSE_LENGTH);
335     switch(*cmd) {
336     case 'D' : // toggle display on/off
337         if (displayOn) {
338             vTaskSuspend(displayHandle);
339             RIT128x96x4Clear();

```

```

340     } else {
341         vTaskResume( displayHandle );
342     }
343     displayOn = !displayOn;
344     ackNack( cData, true );
345
346 #if DEBUG.COMMAND
347     usnprintf( num, 30, "%d %s", displayOn, cData->responseStr );
348     RIT128x96x4StringDraw( num, 0, 30, 15 );
349 #endif
350     break;
351     case 'S' : // start measurements
352         vTaskResume( measureHandle );
353         vTaskResume( ekgCaptureHandle );
354
355         // enable the interrupts used for measurement
356         //     IntEnable( INT_GPIOA ); // for pulse
357         ADCSequenceEnable( ADC0_BASE, 0 ); //     IntEnable( INT_ADC0SS0 ); // for
ekg
358         //     IntEnable( INT_ADC0SS1 ); // for temperature
359         measureOn = true;
360         ackNack( cData, true );
361         break;
362     case 'P' : // stop
363         vTaskSuspend( measureHandle );
364         vTaskSuspend( ekgCaptureHandle );
365
366         // disable the interrupts used for measurement
367         //     IntDisable( INT_GPIOA ); // for pulse
368         ADCSequenceDisable( ADC0_BASE, 0 ); //     IntDisable( INT_ADC0SS0 ); //
for ekg
369         //     IntDisable( INT_ADC0SS1 ); // for temperature
370         measureOn = false;
371         ackNack( cData, true );
372         break;
373     case 'M' : // measure a sensor
374         measureFromSensor( cData );
375         formatResponseStr( cData );
376         break;
377     case 'G' : // Commands for DEBUG mode
378         switch ( *sensor ) {
379             case 'D' :
380                 ackNack( cData, true );
381                 if ( displayOn )
382                     strncat( cData->responseStr, "<p>On</p>", RESPONSE_LENGTH - 1 );
383                 else
384                     strncat( cData->responseStr, "<p>Off </p>", RESPONSE_LENGTH - 1 );
385
386 #if DEBUG.COMMAND
387         usnprintf( num, 30, "%d %s", displayOn, cData->responseStr );
388         RIT128x96x4StringDraw( num, 0, 30, 15 );
389 #endif
390         break;

```



```

391     case 'M' :
392         ackNack(cData, true);
393         addTags((measureOn ? "ON" : "OFF"), false); //TODO fix the false
394     vlaue
395         strncat(cData->responseStr, temporaryBuffer, RESPONSELENGTH - 1);
396 #if DEBUG.COMMAND
397     usnprintf(num, 30, "%d %s", measureOn, cData->responseStr);
398     RIT128x96x4StringDraw(num, 0, 30, 15);
399 #endif
400     break;
401 }
402 ackNack(cData, false);
403 break;
404 default : // send error to remoteStr
405     ackNack(cData, false);
406 #if DEBUG.COMMAND
407     usnprintf(num, 30, "invalid command");
408     RIT128x96x4StringDraw(num, 0, 30, 15);
409 #endif
410 }
411 *(cData->responseReady) = true;
412 vTaskSuspend(NULL);
413 }

```

#### ../code/commandTask.h

```

1  /*
2  * commandTask.h
3  * Author(s): Patrick Ma
4  *
5  * 3/10/2014
6  *
7  * Header for the command task
8  *
9  * Command supported by the commandTask are of the form:
10 * <X> <Y> <Z>
11 *
12 * These placeholders are explained below. If an invalid command sequence is
13 * given, commandTask responds with an 'E' character. Otherwise, the
14 * immediate
15 * response is 'A', followed by any specified return value(s).
16 *
17 * <X> : Primary Command affecting major system functions
18 *
19 * S — Indicates START mode. Causes the hardware to initiate the
20 * previously assigned measurement task.
21 *
22 * P — Indicates STOP mode. Causes the hardware to suspend any current
23 * measurements.
24 *
25 * D — Toggles the local OLED display on and off

```

```

26 *
27 *   M – Initiate a MEASURE function. Takes additional commands to specify
    the
28 *   type and number of measurements
29 *
30 *   G – Causes the commandTask to enter DEBUG mode. Additional modifiers
31 *   specify which state to debug.
32 *
33 *
34 * <Y> : Modifier specifying the sublevel of behavior
35 *
36 *   A – Instruct the system to measure All sensors
37 *
38 *   T – Instruct system to measure Temperature. Will return the
    temperature
39 *   in degrees Celsius
40 *
41 *   B – Instruct system to measure the patient Blood pressure. Will return
42 *   the systolic and diastolic blood pressure in mmHg
43 *
44 *   P – Instruct system to measure patient Pulse rate. Will return the
    heart
45 *   rate in beats per minute
46 *
47 *   E – Instruct the system to perform an EKG measurement. Will return the
48 *   primary EKG frequency
49 *
50 *   S – Query the system battery Status. Returns the percentage left.
51 *
52 *   D – DEBUG mode only. Returns the local Display status
53 *
54 *   M – DEBUG mode only. Returns system Measurement status. Returns
    boolean
55 *   true if measurement is enabled, false otherwise
56 *
57 *   T – DEBUG mode only. Returns the Type of the last measurement
    performed.
58 *
59 *
60 * <Z> : Modifier specifying the number of measurements to perform. The
61 *   measurements are sent to a buffer for usage one at a time.
62 *
63 *   <no value> – Scan continuously
64 *
65 *   <integer value> – perform the specified number of measurements.
66 */
67
68 #include "tcb.h"
69
70 /* points to the TCB for commandTask */
71 extern TCB commandTask;

```

### D.7.12 Status

../code/status.h

```
1  /*
2  * status.h
3  * Author(s): PatrickMa
4  * 1/28/2014
5  *
6  * Defines the public interface for the status task
7  *
8  * initializeStatusTask() should be called once before performing status
9  * functions
10 */
11
12 #include "tcb.h" // for TCBs
13
14 /* Initialize StatusData, must be done before running functions */
15 void initializeStatusTask();
16
17 /* The status Task */
18 extern TCB statusTask;
```

../code/status.c

```
1  /*
2  * status.c
3  * Author(s): PatrickMa
4  * 1/28/2014
5  *
6  * implements status.h
7  */
8
9  #include "status.h"
10 #include "globals.h"
11 #include "timebase.h"
12
13 // StatusData structure internal to compute task
14 typedef struct {
15     unsigned short *batteryState;
16 } StatusData;
17
18 void statusRunFunction(void *data); // prototype for compiler
19
20 static StatusData data; // the internal data
21 TCB statusTask = {&statusRunFunction, &data}; // task interface
22
23 /* Initialize the StatusData task values */
24 void initializeStatusTask() {
25     // Load data
26     data.batteryState = &(global.batteryState);
27
28     // Load TCB
29     statusTask.runTaskFunction = &statusRunFunction;
```

```

30 statusTask.taskDataPtr = &data;
31 }
32
33 /* Perform status tasks */
34 void statusRunFunction(void *data){
35     static tBoolean onFirstRun = true;
36
37     if (onFirstRun) {
38         initializeStatusTask();
39         onFirstRun = false;
40     }
41
42     if (IS_MAJOR_CYCLE) {
43         StatusData *sData = (StatusData *) data;
44         if (*(sData->batteryState) > 0)
45             *(sData->batteryState) = *(sData->batteryState) - 1;
46     }
47 }

```

## D.8 Website Source

../code/fs/index.html

```

1
2 <!DOCTYPE HTML PUBLIC "-//W3C//DTD HTML 4.01 Transitional//EN">
3 <html><head>
4 <meta http-equiv="content-type" content="text/html; charset=ISO-8869-1">
5   title>Doctor at Your Fingertips – Web Interface</title>
6 <script language="JavaScript">
7 <!--
8 function send_command()
9 {
10     var send = false;
11     var command = document.getElementById("terminal");
12
13     // Create a send request (don't bother with ActiveX Objects)
14     send = new XMLHttpRequest();
15
16     if (send)
17     {
18         // User entered non-blank command
19         if (command.value != "")
20         {
21             var res = command.value.split(" ");
22             var uri = encodeURIComponent("/cgi-bin/send_command/value=" + res.join("+"));
23             // Send request to special url for server processing
24             send.open("GET", uri, true);
25             send.send(null);
26             setTimeout(receive_command, 500);
27         }
28     }
29 }

```

```

29 }
30
31 function receive_command() {
32     var receive = false;
33
34     receive = new XMLHttpRequest();
35
36     // Send a request for the current command response
37     if (receive)
38     {
39         receive.open("GET", "/cgi-bin/receive_command", true);
40
41         // Register receive_complete() with this XMLHttpRequest so it is called
42         // when a response is received
43         receive.onreadystatechange = receive_complete;
44         receive.send(null);
45     }
46
47     function receive_complete()
48     {
49         if (receive.readyState == 4)
50         {
51             if (receive.status == 200)
52             {
53                 document.getElementById("response").innerHTML = receive.responseText;
54             }
55         }
56     }
57
58     function blink() {
59         var blinks = document.getElementsByTagName('blink');
60         for (var i = blinks.length - 1; i >= 0; i--) {
61             var s = blinks[i];
62             s.style.visibility = (s.style.visibility === 'visible') ? 'hidden' : 'visible';
63         }
64         window.setTimeout(blink, 1000);
65     }
66     if (document.addEventListener) document.addEventListener("DOMContentLoaded", blink, false);
67     else if (window.addEventListener) window.addEventListener("load", blink, false);
68     else if (window.attachEvent) window.attachEvent("onload", blink);
69     else window.onload = blink;
70     //—>
71 </script>
72 <style type="text/css">
73 body
74 {
75     font-family: Arial;
76     background-color: white;

```

```

77 | margin: 20px;
78 | padding: 0px
79 | }
80 | h1
81 | {
82 |   color: #000000;
83 |   font-family: Arial;
84 |   font-size: 24pt;
85 | }
86 |
87 | </style>
88 | </head>
89 | <body>
90 |
91 | <div id="container" style="float:right; right:50%; position:relative">
92 |   <div id="uicontainer" style="float:right; right:-50%; position:relative">
93 |     <div id="title"> <h1> Doctor at Your Fingertips </h1> </div>
94 |     <div id="ui" style="padding:10px;float:left;border: 1px solid;">
95 |       <div id="enter"> Enter a command: </div>
96 |       <div id="command" style="margin-bottom:10px">
97 |         <input id="terminal" type="text" onkeydown="if (event.keyCode == 13)
          send_command()" />
98 |         <input id="send" value="Send" type="button" onclick="send_command()"
          />
99 |       </div>
100 |     <div id="response"> </div>
101 |   </div>
102 |
103 |   <div id="validcommands" style="margin-left:20%;float:left;border:1px
          solid;">
104 |     <table>
105 |       <tr>
106 |         <td> <b> Command </b> </td>
107 |         <td> <b> Effect/Usage </b> </td>
108 |       </tr>
109 |       <tr>
110 |         <td> S </td>
111 |         <td> Starts the device. Enables measurement devices and
          interrupts. </td>
112 |       </tr>
113 |       <tr>
114 |         <td> P </td>
115 |         <td> Stops the device. Disables measurement devices and
          interrupts. </td>
116 |       </tr>
117 |       <tr>
118 |         <td> D </td>
119 |         <td> Enables/disables the onboard OLED. </td>
120 |       </tr>
121 |       <tr>
122 |         <td> M &lt;sensor(s)&gt; </td>
123 |         <td> Measure the specified sensor(s). Sensors are specified by...
          </td>

```

```
124         </tr>
125     </table>
126 </div>
127 </div>
128 </div>
129
130 </div>
131
132 </body>
133 </html>
```